



HYPEX SURGICAL
INSTRUMENTS

GENERAL INSTRUMENTS CATALOGUE

2026 EDITION

EST. IN SIALKOT, PAKISTAN

**LEADING WORLDWIDE MANUFACTURER
OF SURGICAL & DENTAL INSTRUMENTS**



ISO Aligned



Quality Certified



Batch Traceability



Global Standards

About Us

Hypex Surgical is a Pakistan-based manufacturer of high-quality surgical and dental instruments, built on precision, craftsmanship, and trust. With years of expertise in the medical instruments industry, we take pride in supplying reliable tools that empower healthcare professionals around the world.

Our products are crafted using German-grade stainless steel and advanced production technology, ensuring accuracy, durability, and ergonomic design for every procedure. From general surgery to specialized disciplines, Hypex Surgical instruments reflect our ongoing commitment to excellence and innovation.

We proudly serve global distributors, hospitals, and medical professionals, ensuring timely delivery and unmatched quality standards — because for us, every instrument represents a promise of life-saving precision.

Our **Mission** is to design, manufacture, and deliver world-class surgical instruments that meet the evolving needs of the healthcare industry — blending innovation, craftsmanship, and integrity to support surgeons and improve patient outcomes worldwide.

Hypex Surgical **Vision** is to be recognised globally as a trusted partner in surgical excellence, setting new benchmarks in quality, precision, and reliability. Hypex Surgical aims to continuously evolve with technology, sustainability, and customer satisfaction at its core.

How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work
- 2 Email info@hypexsurgical.com or call +92 332 0640444 with the product codes and we will provide you with a quote
- 3 Upon approval of the quote we will process your order

How To Pay

We accept Visa, Mastercard, Amex cards, paypal and bank payments.



NEED A BESPOKE INSTRUMENT?

CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

Our Services

At Hypex Surgical, we take pride in being a trusted manufacturer and global supplier of high-quality surgical and dental instruments. Each instrument we produce is engineered with precision, built to perform, and designed to last — ensuring reliability in every procedure.

MANUFACTURING EXCELLENCE

We specialize in the design and production of premium-grade surgical and dental instruments using German stainless steel and cutting-edge technology. Our manufacturing process ensures precision, consistency, and durability in every product.

CUSTOM INSTRUMENT DEVELOPMENT

We collaborate closely with surgeons, hospitals, and distributors to develop customized instruments that meet unique procedural needs. From prototyping to final production, every detail is tailored to our client's specifications.

QUALITY CONTROL & ASSURANCE

Every instrument passes through a multi-stage inspection process, including dimensional verification, hardness testing, and surface finishing checks — ensuring compliance with international standards such as ISO 13485 and CE.

GLOBAL EXPORT SERVICES

We proudly serve clients in over 40 countries, providing seamless global logistics and documentation support to ensure smooth and timely deliveries.

PRIVATE LABEL & OEM SERVICES

We provide Private Label and OEM manufacturing options for distributors and brands worldwide, enabling them to market premium instruments under their own identity — supported by our unmatched craftsmanship and reliability.

POLISHING & FINISHING OPTIONS

Choose from a wide range of finishes including matte, mirror, satin, sandblast, and gold-plated, providing both functionality and aesthetic appeal to match your market's demand.

TIMELY DELIVERY & SUPPORT

With a streamlined supply chain and experienced export team, we ensure on-time delivery, transparent communication, and dedicated after-sales support to build long-term trust.

INNOVATION & R&D

Our R&D team continuously explores modern surgical technologies and ergonomic designs to enhance performance, reduce fatigue, and align with the latest trends in global healthcare.

General Instruments


Artery Forceps	04
Clamps	07
Spring Forceps	08
Tissue Forceps	10
Sponge Holding and Sterilising Forceps	12
Other Forceps	13
Scissors	14
Needle Holders	20
Combined Scissor and Needle Holders	23
Bag and Towel Clips	24
Retractors	25
Dissectors and Probes	35
Miscellaneous	36
Silicone Matting	37
Wire Sterilising Baskets	38
Hollow Ware	39



Artery Forceps

Adson Artery Forceps

Box joint, 180mm long



Serrated jaws

Code	Description
HX-GI/001	Straight
HX-GI/001a	Curved




1-2 teeth

Code	Description
HX-GI/001b	Straight
HX-GI/001c	Curved

Cushing Artery Forceps

Box joint, fine, half serrated jaws, 140mm long




Code	Description
HX-GI/002	Straight
HX-GI/002a	Curved



Cairns Artery Forceps

Box joint, half serrated jaws
145mm long



Serrated jaws

Code	Description
HX-GI/003	Straight




1-2 teeth

Code	Description
HX-GI/003a	Curved
HX-GI/003b	Curved to side

Dunhill Artery Forceps

Box joint, fine serrations



Straight


Code	Description
HX-GI/004	135mm long
HX-GI/004a	190mm long

Curved

Code	Description
HX-GI/004b	135mm long
HX-GI/004c	180mm long

Crile Artery Forceps

Box joint, serrated jaws



Straight

Code	Description
HX-GI/005	140mm long
HX-GI/005a	160mm long




Curved

Code	Description
HX-GI/005b	140mm long
HX-GI/005c	160mm long

Lawson Tait Artery Forceps

Box joint, straight jaw




Code	Description
HX-GI/006	140mm long

Artery Forceps

Halstead Mosquito Forceps

Standard pattern, haemostatic, box joint



Serrated jaw

Code	Description
HX-GI/007	Straight 125mm long
HX-GI/007a	Curved 125mm long
HX-GI/007b	Fully curved 125mm long

1-2 teeth

Code	Description
HX-GI/007g	Straight 125mm long
HX-GI/007h	Straight 190mm long
HX-GI/007i	Curved 125mm long
HX-GI/007j	Curved 190mm long

Micro, 100mm long

Code	Description
HX-GI/007c	Straight, serrated
HX-GI/007d	Curved, serrated
HX-GI/007e	Straight, 1-2 teeth
HX-GI/007f	Curved, 1-2 teeth

Kelly Artery Forceps

Box joint, half serrated jaws




Straight

Code	Description
HX-GI/008	140mm long
HX-GI/008a	160mm long

Curved

Code	Description
HX-GI/008b	140mm long
HX-GI/008c	160mm long

Fraser Kelly Artery Forceps Curved




Straight

Code	Description
HX-GI/010	180mm

Heiss Artery Forceps

Box joint, 205mm long



Serrated jaws


Code	Description
HX-GI/009	Straight
HX-GI/009a	Curved

1-2 teeth

Code	Description
HX-GI/009b	Straight
HX-GI/009c	Curved

Kocher Artery Forceps

Box joint, 1-2 teeth



Straight

Code	Description
HX-GI/011	130mm long
HX-GI/011a	150mm long
HX-GI/011b	180mm long
HX-GI/011c	200mm long
HX-GI/011d	230mm long


Curved

Code	Description
HX-GI/011e	150mm long
HX-GI/011f	180mm long
HX-GI/011g	200mm long
HX-GI/011h	230mm long

Artery Forceps

Moynihan Artery Forceps

Box joint, half-serrated jaws



Straight


Code	Description
HX-GI/012	145mm long
HX-GI/012a	180mm long

Curved

Code	Description
HX-GI/012b	145mm long
HX-GI/012c	180mm long

Rochester-Pean Artery Forceps

Box joint



Straight


Code	Description
HX-GI/015	165mm long
HX-GI/015a	190mm long
HX-GI/015b	205mm long
HX-GI/015c	240mm long

Curved

Code	Description
HX-GI/015d	165mm long
HX-GI/015e	180mm long
HX-GI/015f	205mm long
HX-GI/015g	240mm long

Roberts Artery Forceps


Screw joint, curved on flat, serrated



Code	Description
HX-GI/013	230mm long

Spencer Wells Artery Forceps

Box joint, serrated jaws



Straight

Code	Description
HX-GI/016	127mm long
HX-GI/016a	152mm long
HX-GI/016b	177mm long
HX-GI/016c	203mm long


Curved

Code	Description
HX-GI/016d	127mm long
HX-GI/016e	152mm long
HX-GI/016f	177mm long
HX-GI/016g	203mm long

OTHER SIZES AVAILABLE

Dandy Artery Forceps


Box joint, curved in plane



Code	Description
HX-GI/014	140mm long

Carmalt Angiotribe Forceps

Box joint, atraumatic, light pattern



Straight


Code	Description
HX-GI/017	165mm long
HX-GI/017a	200mm long
HX-GI/017b	230mm long

Curved

Code	Description
HX-GI/017c	165mm long

Clamps

Dieffenbach Bulldog Clamp



Straight, serrated

Code	Description
HX-GI/018	1.2mm fine jaw, 40mm long
HX-GI/018a	1.2mm fine jaw, 45mm long
HX-GI/018b	1.8mm jaw, 30mm long
HX-GI/018c	1.8mm jaw, 35mm long
HX-GI/018d	2.5mm jaw, 40mm long
HX-GI/018e	2.5mm jaw, 45mm long
HX-GI/018f	6mm heavy jaw, 60mm long
HX-GI/018g	20mm jaw, 75mm long
HX-GI/018h	30mm jaw, 85mm long
HX-GI/018i	45mm jaw, 105mm long
HX-GI/018j	65mm jaw, 120mm long

Straight, plain

Code	Description
HX-GI/018k	1.2mm fine jaw, 40mm long
HX-GI/018l	1.2mm fine jaw, 45mm long
HX-GI/018m	1.8mm jaw, 30mm long
HX-GI/018n	1.8mm jaw, 35mm long
HX-GI/018o	2.5mm jaw, 40mm long
HX-GI/018p	2.5mm jaw, 45mm long

Curved, serrated

Code	Description
HX-GI/019	1.2mm fine jaw, 40mm long
HX-GI/019a	1.2mm fine jaw, 45mm long
HX-GI/019b	1.8mm jaw, 30mm long
HX-GI/019c	1.8mm jaw, 35mm long
HX-GI/019d	2.5mm jaw, 40mm long
HX-GI/019e	2.5mm jaw, 45mm long
HX-GI/019f	6mm heavy jaw, 60mm long
HX-GI/019g	20mm jaw, 75mm long
HX-GI/019h	30mm jaw, 85mm long
HX-GI/019i	45mm jaw, 105mm long
HX-GI/019j	65mm jaw, 120mm long

Curved, plain

Code	Description
HX-GI/019k	1.2mm fine jaw, 40mm long
HX-GI/019l	1.2mm fine jaw, 45mm long
HX-GI/019m	1.8mm jaw, 30mm long
HX-GI/019n	1.8mm jaw, 35mm long
HX-GI/019o	2.5mm jaw, 40mm long
HX-GI/019p	2.5mm jaw, 45mm long

Blalock Bulldog Clamp



Serrated, with hole in handle

Straight

Code	Description
HX-GI/020	11mm jaw, 38mm long
HX-GI/020a	14mm jaw, 48mm long
HX-GI/020b	17mm jaw, 60mm long
HX-GI/020c	27mm jaw, 70mm long

Curved

Code	Description
HX-GI/020d	11mm jaw, 38mm long
HX-GI/020e	14mm jaw, 48mm long
HX-GI/020f	17mm jaw, 60mm long
HX-GI/020g	27mm jaw, 70mm long

Plain, with hole in handle

Straight

Code	Description
HX-GI/021	11mm jaw, 38mm long
HX-GI/021a	14mm jaw, 48mm long
HX-GI/021b	17mm jaw, 60mm long
HX-GI/021c	27mm jaw, 70mm long

Curved

Code	Description
HX-GI/021d	11mm jaw, 38mm long
HX-GI/021e	14mm jaw, 48mm long
HX-GI/021f	17mm jaw, 60mm long
HX-GI/021g	27mm jaw, 70mm long

An Overview On

Forceps Jaw Types

You will see that forceps in our catalogue will have either a traditional serrated jaw or a different design which is known as atraumatic (or non-traumatic).

Traditional Serrated Jaw



The traditional serrated jaw is generally more robust and used on forceps where the areas are not as delicate.


Atraumatic Jaw



Atraumatic jaws will be mainly found in delicate surgery areas such as vascular and cardi thoracic. The design is such, that the surface area gripped is reduced, ensuring less trauma without the loss of grip.

Spring Forceps


Adsons Dissecting Forceps



127mm long

Code	Description
HX-GI/022	1-2 teeth
HX-GI/022a	Lightweight, micro
HX-GI/022b	Serrated
HX-GI/022c	Serrated, lightweight, micro


Block End Dissecting Forceps



Serrated

Code	Description
HX-GI/023	127mm long
HX-GI/023a	152mm long
HX-GI/023b	203mm long

Bonney Dissecting Forceps




178mm long

Code	Description
SC-GI/267	Serrated, 1-2 teeth
SC-GI/268	1-2 teeth

Charnley-McIndoe Dissecting Forceps


150mm long



Code	Description
HX-GI/025	2-3 teeth

DeBakey Atraumatic Dissecting Forceps

Straight, atraumatic jaws



1.5mm

Code	Description
HX-GI/026	150mm long
HX-GI/026a	240mm long

2.5mm

Code	Description
HX-GI/026d	150mm long
HX-GI/026e	240mm long

2mm


Code	Description
HX-GI/026b	150mm long
HX-GI/026c	180mm long

2.8mm

Code	Description
HX-GI/026f	150mm long
HX-GI/026g	240mm long
HX-GI/026h	180mm long

Dissecting Forceps

Turnover end, serrated jaws



Code	Description
HX-GI/027	127mm long
HX-GI/027a	152mm long

Code	Description
HX-GI/027b	178mm long
HX-GI/027c	203mm long

Dunhill Dissecting Forceps


Light, fine serrated jaws



Code	Description
HX-GI/028	145mm long

Fredrick Coole Dissecting Forceps

Extra fine serrated jaws




125mm long

Code	Description
HX-GI/029	Straight
HX-GI/029a	Curved

Spring Forceps

Gillies Dissecting Forceps


1-2 teeth



Code	Description
HX-GI/030	152mm long

McIndoe Dissecting Forceps


152mm long



Code	Description
HX-GI/034	Serrated
HX-GI/034a	1-2 teeth

Jackson Burrows Dissecting Forceps


140mm long



Code	Description
HX-GI/031	Tongue and groove, 1-2 teeth
HX-GI/031a	Tongue and groove only

McIndoe Modified Dissecting Forceps


152mm long, 1mm tip, modified for fine powerful dissection with strong action tip



Code	Description
HX-GI/035	Serrated


Jean Dissecting Forceps

Serrated jaws, 140mm long



Code	Description
HX-GI/032	Fine
HX-GI/032a	1-2 teeth
HX-GI/032b	2-3 teeth, 178mm long

Treves Dissecting Forceps



1-2 teeth


Code	Description
HX-GI/036	127mm long
HX-GI/036a	152mm long

2-3 teeth

Code	Description
HX-GI/036b	127mm long
HX-GI/036c	152mm long

Lanes Dissecting Forceps


1-2 teeth on serrated jaws, 140mm long



Code	Description
HX-GI/033	127mm long
HX-GI/033a	152mm long
HX-GI/033b	178mm long

Code	Description
HX-GI/033c	2-3 teeth, 178mm long
HX-GI/033d	2-3 teeth 203mm long

Waugh Dissecting Forceps



Fine, serrated

Code	Description
HX-GI/037	152mm long
HX-GI/037a	178mm long
HX-GI/037b	203mm long


Fine, 1-2 teeth

Code	Description
HX-GI/037c	152mm long
HX-GI/037d	178mm long
HX-GI/037e	203mm long

Tissue Forceps

Allis Tissue Forceps

Box joint



150mm long


Code	Description
HX-GI/038	3-4 teeth
HX-GI/038a	4-5 teeth

190mm long

Code	Description
HX-GI/038b	5-6 teeth

Duval Tissue Forceps

Box joint




Code	Description
HX-GI/041	13mm wide, 203mm long

Littlewood Tissue Forceps

Box joint, 2-3 teeth


Code	Description
HX-GI/039	190mm long



Lane Tissue Forceps

Box joint, 1-2 teeth, fenestrated


Code	Description
HX-GI/042	127mm long
HX-GI/042a	152mm long
HX-GI/042b	197mm long



Babcock Tissue Forceps

Box joint, serrated jaws


Code	Description
HX-GI/040	160mm long
HX-GI/040a	200mm long
HX-GI/040b	230mm long



Rutherford Morrison Tissue Forceps

Box joint


Code	Description
HX-GI/043	4-5 teeth, 152mm long



Tissue Forceps

Shaw Tissue Forceps

Box joint




Code	Description
HX-GI/044	1-2 teeth, 178mm long



Guys Gallabin Suture-Holding Forceps


Screw joint, angled flat blades



Code	Description
HX-GI/045	125mm long
HX-GI/045a	150mm long

Russian Dissecting Forceps

With fenestrated end



Code	Description
HX-GI/046	150mm long
HX-GI/046a	200mm long
HX-GI/046b	230mm long
HX-GI/046c	250mm long




Special Offer

on SSD Essentials Pack


Top up your trays with these essentials...

CONTACT
US FOR
OUR OFFER
PRICE

Halstead Mosquito Forceps
Straight/Curved
125mm (5")




Rampléy Sponge Holding Forceps
Straight, Serrated
Jaws 241mm
(9½")




Scalpel Blade Handle
No. 3



Spring Clip
Serrated Pattern,
Stainless Steel for
Fixing Slings to Leg
Splints, also for
Towel Clipping
6.3cm (2½") Wide



Spencer Wells Artery Forceps
Straight, Serrated
Jaws




Email us to place your order

info@hypexsurgical.com

Sponge Holding and Sterilising Forceps

Rampley Sponge Holding Forceps


Box joint, straight, serrated jaws




Code	Description
HX-GI/047	178mm long
HX-GI/047a	241mm long

Lister Sinus Forceps

Box joint, serrated tip




Code	Description
HX-GI/050	127mm long
HX-GI/050a	152mm long
HX-GI/050b	178mm long



Foerster Sponge Holding Forceps

Box joint, straight, serrated jaws




Code	Description
HX-GI/048	178mm long



Cheatle Sterilizer Forceps

Screw joint




Code	Description
HX-GI/051	267mm long

Harrison Bowl Steriliser Forceps

Screw joint, angled, double jaw

Code	Description
HX-GI/049	356mm long
HX-GI/049a	457mm long



Other Forceps

Hunter Splinter Forceps


Serrated grooved pointed jaws



Code	Description
HX-GI/052	115mm long

Michel Forceps


Box joint, double-ended, for inserting and removing clips stainless steel



Code	Description
HX-GI/055	130mm long

Michel Skin Holding Forceps


With magazine for Michel clips



Code	Description
HX-GI/053	1-2 teeth, 127mm long

Michel Clip Extracting Forceps


Box joint



Code	Description
HX-GI/056	125mm long

Childe Approximating Forceps


With magazine for Michel clips



Code	Description
HX-GI/054	180mm long

Robins Anchoring Forceps


140mm long



Code	Description
HX-GI/057	5mm with holder for diathermy leads
HX-GI/057a	13mm with holder for suction tubing

Premium Quality Stainless Steel Scissors

Dressing Scissors



Blunt/Blunt Straight

Code	Description
HX-GI/058	130mm long
HX-GI/058a	146mm long
HX-GI/058b	165mm long

Curved

Code	Description
HX-GI/058f	130mm long
HX-GI/058g	146mm long
HX-GI/058h	165mm long
HX-GI/058i	185mm long
HX-GI/058j	147mm long

With one serrated blade

Blunt/Sharp Straight

Code	Description
HX-GI/058c	130mm long
HX-GI/058d	146mm long
HX-GI/058e	165mm long

Sharp/Sharp Straight

Code	Description
HX-GI/058k	130mm long
HX-GI/058l	146mm long
HX-GI/058m	165mm long

Curved


Code	Description
HX-GI/058n	130mm long

Dressing Scissors, Left-Handed

With one black bow, 145mm long

Blunt/Blunt

Code	Description
HX-GI/059	Straight
HX-GI/059a	Curved



Blunt/Sharp


Code	Description
HX-GI/059b	Straight
HX-GI/059c	Curved

Sharp/Sharp

Code	Description
HX-GI/059d	Straight
HX-GI/059e	Curved

Littler Dissecting Scissors

Fine pattern, with eye for suture



Code	Description
HX-GI/060	115mm long


Kilner Dissecting Scissors

Fine pattern



Code	Description
HX-GI/061	115mm long

Metzenbaum Scissors



Straight

Code	Description
HX-GI/062	180mm long
HX-GI/062a	147mm long
HX-GI/062b	200mm long
HX-GI/062c	230mm long
HX-GI/062d	250mm long
HX-GI/062e	280mm long

Curved

Code	Description
HX-GI/062f	185mm long
HX-GI/062g	147mm long
HX-GI/062h	200mm long
HX-GI/062i	230mm long
HX-GI/062j	250mm long
HX-GI/062k	280mm long

Mayo Scissors



Straight


Code	Description
HX-GI/063	140mm long
HX-GI/063a	160mm long
HX-GI/063b	190mm long
HX-GI/063c	230mm long

Curved

Code	Description
HX-GI/063d	140mm long
HX-GI/063e	160mm long
HX-GI/063f	190mm long
HX-GI/063g	230mm long


Fergusson Scissors

Angled on flat



Code	Description
HX-GI/064	180mm long


Paynes Suture Scissors



Code	Description
HX-GI/065	100mm long

Ligature Stitch Scissors

Fine pattern




Code	Description
HX-GI/066	115mm long



FOR INFORMATION ON
BLACK AND GOLD BOWS
SEE PAGE 17

Scissors


Lister Bandage Scissors



Code	Description
HX-GI/067	147mm long
HX-GI/067a	178mm long

Lister Bandage Scissors


Angled, with probe end



Code	Description
HX-GI/071	147mm long
HX-GI/071a	178mm long

McIndoe Scissors

178mm long




Code	Description
HX-GI/068	Straight
HX-GI/068a	Curved on flat
HX-GI/068b	Supercut, curved




Thomas Strabismus Scissors

114mm long




Code	Description
HX-GI/072	Straight
HX-GI/072a	Curved

Heath Stitch Scissors




Code	Description
HX-GI/069	152mm long

"Careless" Scissors



Code	Description
HX-GI/073	165mm long

Tubing Scissors



Code	Description
HX-GI/070	130mm long

Utility Scissors




Code	Description
HX-GI/074	190mm long

Tungsten Carbide Hard Insert Scissors

Dressing Scissors

Tungsten carbide hard insert, blunt/blunt



Straight


Code	Description
HX-GI/075	130mm long
HX-GI/075a	145mm long
HX-GI/075b	165mm long
HX-GI/075c	185mm long

Curved

Code	Description
HX-GI/075d	130mm long
HX-GI/075e	145mm long
HX-GI/075f	165mm long

Mayo Dressing Scissors

Tungsten carbide hard insert



Straight


Code	Description
HX-GI/079	140mm long
HX-GI/079a	160mm long
HX-GI/079b	190mm long
HX-GI/079c	230mm long

Curved

Code	Description
HX-GI/079d	140mm long
HX-GI/079e	160mm long
HX-GI/079f	190mm long
HX-GI/079g	230mm long

Dressing Scissors

Tungsten carbide hard insert, blunt/sharp



Straight


Code	Description
HX-GI/076	130mm long
HX-GI/076a	145mm long
HX-GI/076b	165mm long
HX-GI/076c	185mm long

Curved

Code	Description
HX-GI/076d	130mm long
HX-GI/076e	145mm long
HX-GI/076f	165mm long

Mayo-Harrington Dressing Scissors

Tungsten carbide hard insert



230mm long


Code	Description
HX-GI/080	Straight
HX-GI/080a	Curved

200mm long

Code	Description
HX-GI/080b	Straight
HX-GI/080c	Curved

Dressing Scissors

Tungsten carbide hard insert sharp/sharp



Straight


Code	Description
HX-GI/077	130mm long
HX-GI/077a	145mm long
HX-GI/077b	165mm long
HX-GI/077c	185mm long

Curved

Code	Description
HX-GI/077d	130mm long
HX-GI/077e	145mm long
HX-GI/077f	165mm long

Metzenbaum Scissors

Tungsten carbide hard insert



Curved


Code	Description
HX-GI/081	145mm long
HX-GI/081a	185mm long
HX-GI/081b	230mm long
HX-GI/081c	250mm long
HX-GI/081d	280mm long

Straight

Code	Description
HX-GI/081e	145mm long
HX-GI/081f	185mm long
HX-GI/081g	200mm long
HX-GI/081h	230mm long
HX-GI/081i	250mm long
HX-GI/081j	280mm long

Baby Metzenbaum Scissors


Tungsten carbide, hard insert, 145mm long



Code	Description
HX-GI/078	Straight
HX-GI/078a	Curved

Ligature Scissors

Tungsten carbide hard insert, curved




Code	Description
HX-GI/082	145mm long
HX-GI/082a	185mm long

Tungsten Carbide Supercut Scissors

Mayo Harrington Scissors

Supercut, tungsten carbide



230mm long


Code	Description
HX-GI/083	Straight
HX-GI/083a	Curved

200mm long

Code	Description
HX-GI/083b	Straight
HX-GI/083c	Curved

Mayo Scissors

Supercut, tungsten carbide



Straight

Code	Description
HX-GI/084	140mm long
HX-GI/084a	160mm long
HX-GI/084b	190mm long
HX-GI/084c	230mm long

Curved

Code	Description
HX-GI/084d	140mm long
HX-GI/084e	160mm long
HX-GI/084f	190mm long
HX-GI/084g	230mm long

Metzenbaum Dissecting Scissors

Supercut, tungsten carbide



Straight

Code	Description
HX-GI/085	185mm long
HX-GI/085a	200mm long
HX-GI/085b	230mm long
HX-GI/085c	250mm long
HX-GI/085d	280mm long

Curved

Code	Description
HX-GI/085e	147mm
HX-GI/085f	185mm
HX-GI/085g	200mm
HX-GI/085h	230mm
HX-GI/085i	250mm
HX-GI/085j	280mm

Dressing Scissors

Supercut, tungsten carbide



Blunt/Blunt

Straight

Code	Description
HX-GI/086	130mm long
HX-GI/086a	145mm long
HX-GI/086b	165mm long
HX-GI/086c	185mm long

Curved

Code	Description
HX-GI/086d	130mm long
HX-GI/086e	145mm long
HX-GI/086f	165mm long
HX-GI/086g	185mm long

Blunt/Sharp

Straight

Code	Description
HX-GI/086h	130mm long
HX-GI/086i	145mm long
HX-GI/086j	165mm long
HX-GI/086k	185mm long

Curved

Code	Description
HX-GI/086l	130mm long
HX-GI/086m	145mm long
HX-GI/086n	165mm long
HX-GI/086o	185mm long

Sharp/Sharp

Straight

Code	Description
HX-GI/086p	130mm long
HX-GI/086q	145mm long
HX-GI/086r	165mm long
HX-GI/086s	185mm long

Curved

Code	Description
HX-GI/086t	130mm long
HX-GI/086u	145mm long
HX-GI/086v	165mm long
HX-GI/086w	185mm long

Guide for Gold and Black Bows

It's important to know the difference for the intended use of each scissor:



2 Gold Bows
Tungsten carbide insert
or padding for a scissor



2 Black Bows
Supercut edge



1 Black and 1 Gold Bow
Supercut tungsten
carbide



1 Black Bow (left ring)
Left-handed instrument




WITH TWO BLACK BOWS
FOR IDENTIFICATION

Supercut Operating Scissors

Dressing Scissors

Supercut, blunt/blunt




145mm long

Code	Description
HX-GI/087	Straight
HX-GI/087a	Curved

Dressing Scissors

Supercut, blunt/sharp




145mm long

Code	Description
HX-GI/088	Straight
HX-GI/088a	Curved

Dressing Scissors

Supercut, sharp/sharp




145mm long

Code	Description
HX-GI/089	Straight
HX-GI/089a	Curved

Mayo Scissors

Supercut



Straight


Code	Description
HX-GI/090	140mm long
HX-GI/090a	160mm long
HX-GI/090b	170mm long

Curved

Code	Description
HX-GI/090c	140mm long
HX-GI/090d	160mm long
HX-GI/090e	170mm long

Boyd Scissors


Supercut, blunt/blunt



Code	Description
HX-GI/091	180mm long

Metzenbaum Dissecting Scissors

Supercut



Straight (length)

Code	Description
HX-GI/092	147mm
HX-GI/092a	185mm
HX-GI/092b	200mm
HX-GI/092c	230mm
HX-GI/092d	250mm
HX-GI/092e	280mm

Curved (length)

Code	Description
HX-GI/092f	147mm
HX-GI/092g	185mm
HX-GI/092h	200mm
HX-GI/092i	230mm
HX-GI/092j	250mm
HX-GI/092k	280mm

Stevens Tenotomy Scissors

Supercut



Curved/blunt

Code	Description
HX-GI/093	114mm long
HX-GI/093a	125mm long
HX-GI/093b	135mm long

Curved/sharp

Code	Description
HX-GI/093c	114mm long
HX-GI/093d	140mm long
HX-GI/093e	160mm long


Code	Description
HX-GI/093f	Straight/blunt, 114mm long
HX-GI/093g	Straight/sharp, 114mm long
HX-GI/093h	Straight/sharp, 180mm long

With cutout, 110mm long

Code	Description
HX-GI/093i	Straight/blunt
HX-GI/093j	Curved/blunt

Strabismus Scissors

Supercut




114mm long

Code	Description
HX-GI/094	Straight
HX-GI/094a	Curved

Iris Scissors

Supercut




114mm long

Code	Description
HX-GI/095	Straight
HX-GI/095a	Curved

Needle Holders

Bruce Clarke Needle Holder


Screw joint



Code	Description
HX-GI/096	127mm long

Lawrence Needle Holder


Fine jaw, 150mm long



Code	Description
HX-GI/099	Fine jaw, 150mm long
HX-GI/099a	Tungsten carbide jaw

Crile-Wood Needle Holder


Straight, baby, fine



Code	Description
HX-GI/097	150mm long

Needle Holder


Surgical Holdings pattern, box joint, lightweight, with curved plain jaws and one offset bow for thumb



Code	Description
HX-GI/100	152mm long

Kilner Needle Holder


Box joint, cross-serrated jaw



Code	Description
HX-GI/098	165mm long
HX-GI/098a	133mm long

MacPhail Needle Holder

Screw joint, with one copper jaw




Code	Description
HX-GI/101	159mm long
HX-GI/101a	178mm long
HX-GI/101b	209mm long

Needle Holders

Mathieu Needle Holder


Screw joint



Code	Description
HX-GI/102	172mm long
HX-GI/102a	209mm long
HX-GI/102b	243mm long



Naughton Morgan Needle Holder

Box joints, cross serrated jaws



Mayo Needle Holder

Box joint, straight



Wide cross serrated jaws

Code	Description
HX-GI/103	159mm long
HX-GI/103a	203mm long

Narrow cross serrated jaws

Code	Description
HX-GI/103b	159mm long


Tungsten carbide

Code	Description
HX-GI/104	185mm long
HX-GI/104a	254mm long

Code	Description
HX-GI/104b	185mm long
HX-GI/104c	254mm long

Rochester-Mayo Needle Holder

Screw joint



Code	Description
HX-GI/105	178mm long

Choosing the most suitable needle size

Our tungsten carbide inserts are available in a variety of styles. Choosing the most suitable needle size for each procedure helps give the surgeon maximum grip and control.

Coarse 0.65mm

For powerful needle holder patterns and wire twisters. Maximum holding capabilities for heavy needles or wire.



Standard 0.5mm

Used, for example, in Mayo-Hegar, Mathieu Needle Holders. For needles 3-0 and larger.



Fine 0.4mm

For delicate needle holders such as Crile-Wood and DeBakey. For needles and sutures between 4-0 and 6-0.



Micro 0.2-0.3mm

For micro needle holders such as Castroviejo, Microvascular and Mini-Ryder. For delicate needles and sutures between 7-0 and 10-0.



Smooth

Smooth tungsten carbide surface reduces damage to delicate needles and suture of 6-0 and finer.






WITH TWO GOLD BOWS
FOR IDENTIFICATION

Tungsten Carbide Needle Holders

Berry Needle Holder/Wire Twister


Tungsten carbide jaws



Code	Description
HX-GI/106	200mm long

Finochietto Needle Holder


Tungsten carbide jaws



Code	Description
HX-GI/109	215mm long
HX-GI/109a	265mm long

Crile-Wood Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



Code	Description
HX-GI/107	159mm long
HX-GI/107a	178mm long
HX-GI/107b	150mm long

Halsey Needle Holder

Tungsten carbide insert, box joint, satin finish, straight



Code	Description
HX-GI/110	127mm long

Derf Needle Holder


Tungsten carbide insert



Code	Description
HX-GI/108	115mm long

Kilner Needle Holder

Tungsten carbide insert




Code	Description
HX-GI/111	140mm long



Tungsten Carbide Needle Holders

Masing Needle Holder


Tungsten carbide insert,
one angled bow



Code	Description
HX-GI/112	130mm long

Nievert Needle Holder


Tungsten carbide insert, box joint,
satin finish, straight



Code	Description
HX-GI/113	127mm long

DeBaKey Needle Holder


250mm long



Code
HX-GI/114

Ryder Needle Holder

Tungsten carbide insert, box joint, satin finish, straight



0.2mm jaw pitch


Code	Description
HX-GI/115	135mm long
HX-GI/115a	155mm long
HX-GI/115b	175mm long
HX-GI/115c	210mm long

0.4mm jaw pitch

Code	Description
HX-GI/115d	135mm long
HX-GI/115e	155mm long
HX-GI/115f	175mm long
HX-GI/115g	210mm long

Sarot Needle Holder

Diamond jaw, tungsten carbide, gold plated bows



Code	Description
HX-GI/116	181mm long

Tungsten Carbide Needle Holders

Features and Benefits

Jaw

Tungsten carbide inserts
bonded into the jaws for
a long lasting, secure grip

Bevelled Edges

Prevents suture from snagging

Laser Etching

Date and catalogue number
are etched for easy re-ordering.
Laser etch will not harm
instrument like traditional
methods

Satin Finish

Reduces glare

Ratchets

Precisely cut for effortless
locking/unlocking

Teeth

Fine, hard teeth, diamond
cut into an insert of hardened
tungsten carbide for superior
needle grasping and suture
tying capabilities

Vacuum Brazing

Combines bonding of the
tungsten carbide insert with
the simultaneous vacuum
hardening of the instrument
for greater durability and
longevity

Shanks


Precision balanced for
comfort and control

Gold Handles

For easy identification

Mayo-Hegar Needle Holder

Tungsten carbide insert, box joint, satin finish, straight




Code	Description
HX-GI/117	130mm long
HX-GI/117a	159mm long
HX-GI/117b	178mm long
HX-GI/117c	209mm long

LONGER SIZES AVAILABLE

Combined Scissor and Needle Holders

Olsen-Hegar Scissor/Needle Holder


With tungsten carbide insert



Code	Description
HX-GI/118	147mm long
HX-GI/118a	165mm long

Olsen-Hegar Scissor/Needle Holder


Screw joint



Code	Description
HX-GI/121	140mm long
HX-GI/121a	165mm long
HX-GI/121b	190mm long

Gillies Scissor/Needle Holder


Right-hand, with tungsten carbide insert



Code	Description
HX-GI/119	158mm long

Gillies Scissor/Needle Holder


Heavy pattern, 158mm long



Code	Description
HX-GI/122	Right-hand
HX-GI/122a	Left-hand

Gillies Scissor/Needle Holder


158mm long



Code	Description
HX-GI/120	Right-hand
HX-GI/120a	Left-hand

Foster Gillies Scissor/Needle Holder

Fine pattern




Code	Description
HX-GI/123	127mm long

Bag and Towel Clips

Spring Clip

Serrated pattern, stainless steel, for fixing slings to leg splints, also for towel clipping



Code	Description
HX-GI/124	63mm wide

Spring Bag Clip


Serrated



Code	Description
HX-GI/125	26mm wide

Jones Towel Clip


Superior quality spring forcep type



Code	Description
HX-GI/126	89mm long

Shardles Towel Clip


Spring curved, cross action



Code	Description
HX-GI/127	77mm long

Backhaus Towel Clip


Superior quality, box joint



Code	Description
HX-GI/128	89mm long
HX-GI/128a	102mm long
HX-GI/128b	130mm long

Moynihan Towel Clip


Screw joint, 2-2 teeth, for sides of incision



Code	Description
HX-GI/129	191mm long

Ball and Socket Towel Forceps

Box joint




Code	Description
HX-GI/130	89mm long
HX-GI/130a	115mm long
HX-GI/130b	140mm long

Retractors

Goligher Retractor

Complete with 2 centre blades:
1 x 90mm wide x 50mm deep
1 x 90mm wide x 90mm deep

and 2 pairs of lateral blades:
1 pair 25mm wide x 45mm deep
1 pair 50mm wide x 90mm deep



Code

HX-GI/131

Balfour Abdominal Retractor

Trivalve, frame 203mm

Child


102mm opening, 51mm deep lateral blades,
38mm deep centre blade

Adult

178mm opening, 64mm deep lateral blades,
44mm deep centre blade

Extra large

178mm opening, 102mm deep lateral blades,
108mm deep centre blade




Code	Description
HX-GI/132	Child
HX-GI/132a	Adult
HX-GI/132b	Extra large

Miles Retractor

Complete with 2 centre blades:
1 x 90mm wide x 50mm deep
1 x 90mm wide x 75mm deep

and 2 pairs of lateral blades:
1 pair 80mm wide x 35mm deep
1 pair 80mm wide x 75mm deep




Code	Description
HX-GI/133	340mm long

OTHER BLADE SIZES AVAILABLE

Balfour Abdominal Retractor


Bivalve, baby



Code	Description
HX-GI/134	Fully open 89mm

Berkeley Abdominal Retractor

Parallel action




Code	Description
HX-GI/135	Blades 35mm x 60mm deep, 215mm long

Retractors

Collin Retractor


110mm fully open, 235mm long



Code	Description
HX-GI/136	Complete with centre blade 55 x 40mm and one pair of blades 38mm x 60mm
HX-GI/136a	Without centre blade

Gelpi Retractor


Self retaining



Code	Description
HX-GI/137	140mm long
HX-GI/137a	178mm long


Joll Thyroid Retractor

Self retaining, with lock-in nut



Code	Description
HX-GI/138	152mm long


Gossett Abdominal Retractor




Code	Description
HX-GI/139	Large, fully open 185mm long
HX-GI/139a	Small, fully open 95mm long

Retractors

Pozzi Abdominal Retractors




Code	Description
HX-GI/140	Swivel blades: 38 x 60mm, bivalve, 254mm long



Code	Description
HX-GI/140a	Swivel blades: 38 x 59mm, trivalve, 203mm long

Bowlby Dilator


Curved



Code	Description
HX-GI/141	145mm long

Travers Retractor

Self retaining with internal rack, 210mm long




Code	Description
HX-GI/142	4 x 5 teeth, 95mm fully open

Retractors

Irwin Retractor


Self retaining, with cam rack, 4 x 4 teeth



Code	Description
HX-GI/143	178mm long

Norfolk and Norwich Retractor


4 x 5 teeth 50mm blades x 223mm long



Code	Description
HX-GI/144	Standard rack

West Retractor


Self retaining, 140mm long



Code	Description
HX-GI/145	3 x 4 teeth, 60mm fully open
HX-GI/145a	3 x 4 teeth, 60mm fully open, with cam rack

Weitlander Retractor

3-4 teeth



135mm long


Code	Description
HX-GI/146	Blunt
HX-GI/146a	Sharp
HX-GI/146b	Sharp, 2-3 teeth
HX-GI/146c	Blunt, 2-3 teeth

165mm long

Code	Description
HX-GI/146d	Blunt
HX-GI/146e	Sharp

Retractors

Kocher Retractor



230mm long


Code	Description
HX-GI/148	36 x 20mm blade

250mm long

Code	Description
HX-GI/148a	67 x 25mm blade
HX-GI/148b	80 x 42mm blade
HX-GI/148c	80 x 50mm blade
HX-GI/148d	80 x 60mm blade

Kocher Langenbeck Retractor

215mm long




Code	Description
HX-GI/149	25 x 6mm blade
HX-GI/149a	35 x 8mm blade
HX-GI/149b	35 x 11mm blade
HX-GI/149c	35 x 15mm blade
HX-GI/149d	41 x 11mm blade

Code	Description
HX-GI/149e	55 x 11mm blade
HX-GI/149f	70 x 15mm blade
HX-GI/149g	80 x 12mm blade
HX-GI/149h	80 x 16mm blade


Langenbeck Retractor

Single ended



Code	Description
HX-GI/150	Blade 22mm x 6mm
HX-GI/150a	Blade 45mm x 13mm
HX-GI/150b	Blade 65mm x 13mm
HX-GI/150c	Blade 75mm x 25mm
HX-GI/150d	Blade 65mm x 25mm

Landon Retractor




Code	Description
HX-GI/151	230mm



Canny Ryall Navy Retractor


Double-ended, with one single hook blade 5mm wide, and one fenestrated blade



Code	Description
HX-GI/152	Fenestrated blade 16mm wide, 190mm long
HX-GI/152a	Fenestrated blade 19mm wide, 190mm long
HX-GI/152b	Fenestrated blade 25mm wide, 220mm long
HX-GI/152c	Fenestrated blade 30mm wide, 220mm long

Canny Ryall Army Retractor

Double prong, 13mm wide, 165mm long




Code	Description
HX-GI/153	Large, 19mm wide
HX-GI/153a	Small, 11mm wide



Retractors


Czerny Retractor

Double ended



Code	Description
HX-GI/154	171mm long


Elliot (Modified Deaver) Retractor



Code	Description
HX-GI/157	With flat handle


Lane Retractor

Double ended blade, 229mm long




Code	Description
HX-GI/155	13mm x 19mm wide
HX-GI/155a	25mm x 32mm wide

Deaver Abdominal Retractors



Code	Description
HX-GI/158	Set of 5, 315mm overall length
HX-GI/158a	25mm wide
HX-GI/158b	38mm wide
HX-GI/158c	51mm wide
HX-GI/158d	76mm wide


Moynihan Retractor



Code	Description
HX-GI/156	Single ended, blade 13 x 22mm, 215mm long

Deaver Abdominal Retractors

With round handles




Code	Description
HX-GI/159	Set of 5, 305mm overall length
HX-GI/159a	25mm wide
HX-GI/159b	25mm wide
HX-GI/159c	38mm wide
HX-GI/159d	51mm wide
HX-GI/159e	76mm wide

Retractors

Morris Retractor

Single blade, 229mm long



Code	Description
SC-GI/330	38 x 51mm
SC-GI/331	51 x 51mm
SC-GI/332	64 x 51mm

Kelly Retractor Blade


267mm long



Code	Description
HX-GI/163	38 x 191mm
HX-GI/163a	57 x 191mm

Morris Retractor


Double ended, 216mm long



Code	Description
SC-GI/327	32 x 51mm and 38 x 51mm
SC-GI/328	44 x 51mm and 51 x 51mm
SC-GI/329	64 x 51mm and 38 x 51mm

Durham Retractor


With fenestrated blade, 216mm long



Code	Description
HX-GI/164	10mm wide
HX-GI/164a	14mm wide
HX-GI/164b	20mm wide

Dyball and Lock Retractor

For liver



Code	Description
HX-GI/162	64 x 144mm, 241mm long

OTHER SIZES AVAILABLE

Ballance (Volkman) Retractor

203mm long



Sharp


Code	Description
HX-GI/165	1 prong
HX-GI/165a	2 prong
HX-GI/165b	3 prong
HX-GI/165c	4 prong
HX-GI/165d	6 prong

Blunt

Code	Description
HX-GI/165e	1 prong
HX-GI/165f	2 prong
HX-GI/165g	3 prong
HX-GI/165h	4 prong
HX-GI/165i	6 prong

Retractors

Volkman Retractor



Code	Description
HX-GI/166	Sharp, 220mm long
HX-GI/166a	Blunt, 220mm long
HX-GI/166b	Sharp, 245mm long

Kilner Retractor

Double ended, cat's paw



Code	Description
HX-GI/167	152mm long

Kocher Retractor



Code	Description
SC-OR/N036c	Sharp, 205mm long
SC-OR/N036	Sharp, 220mm long
SC-OR/N036b	Blunt, 220mm long

Senn-Miller Retractor

170mm long



Code
HX-GI/169

Single Hook Retractor

Sharp



Code	Description
HX-GI/170	152mm long
HX-GI/170a	Half size, fine

Gillies Skin Hook

Solid handle, 190mm long




Code	Description
HX-GI/171	Standard, fine
HX-GI/171a	Half size, fine

Turner Warwick Retractors

Turner Warwick Universal Self Retaining Retractor

Code	Description
HX-GI/172	With frames A/B/C/D and blades A/B/C/D/E/F/K/L
HX-GI/172a	Circular frame A only 178mm
HX-GI/172b	Circular frame B only 229mm
HX-GI/172c	Circular frame C only 305mm
HX-GI/172d	Oval frame D only 475mm x 330mm
HX-GI/172e	Oval frame E only 370mm x 300mm



Turner Warwick Needleholder with Tungsten Carbide Jaws

Curved shanks with one short bow


Code	Description
HX-GI/173	230mm long
HX-GI/173a	200mm long

Blade A




Code	Description
HX-GI/174a	Blade only - 19mm x 76mm

Blade B




Code	Description
HX-GI/174b	Blade only - 38mm x 102mm

Blade C



Code	Description
HX-GI/174c	Blade only - 57mm x 89mm


Blade D



Code	Description
HX-GI/174d	Blade only - 57mm x 114mm


Turner Warwick Retractors

Blade E




Code	Description
HX-GI/174e	Blade only - curved stainless steel blade 51mm x 38mm

Blade F



Code	Description
HX-GI/174f	Blade only - curved stainless steel blade 76mm x 60mm

Blade G




Code	Description
HX-GI/174g	Blade only - flat stainless steel blade 89mm x 60mm

Blade J




Code	Description
HX-GI/174j	Blade only - single hook, 13mm deep

Blade K



Code	Description
HX-GI/174k	Blade only - flat with clip 55mm x 90mm

Blade L



Code	Description
HX-GI/174l	Blade only - flat with clip 55mm x 120mm

Dissectors and Probes

MacDonald Dissector Double Ended



Code	Description
HX-GI/175	191mm long

Probe, Sinus Silver



Code	Description
HX-GI/179	152mm long

Brodie Director Probe Ended, Silver Plated



Code	Description
HX-GI/176	165mm long

Probe, Fish Tail Silver Plated



Code	Description
HX-GI/180	127mm long
HX-GI/180a	152mm long
HX-GI/180b	203mm long
HX-GI/180c	254mm long

Watson Cheyne Dissector

Fine, double-ended, with probe end



Code	Description
HX-GI/177	127mm long
HX-GI/177a	191mm long
HX-GI/177b	230mm long
HX-GI/177c	255mm long

Jobson Horne Probe Dressing (Aural)

With one flat tapered end, serrated edge and one small ring end, stainless steel



Code	Description
HX-GI/181	140mm long
HX-GI/181a	178mm long

Probe, Malleable Silver Plated



Code	Description
HX-GI/178	127mm long
HX-GI/178a	152mm long
HX-GI/178b	178mm long
HX-GI/178c	203mm long
HX-GI/178d	255mm long
HX-GI/178e	305mm long

Miscellaneous


Aneurysm Needle

Small, for infusions




Code	Description
HX-GI/182	146mm long

Syme Aneurysm Needle



Code	Description
HX-GI/183	165mm long


Ring Cutter



Code	Description
HX-GI/184	Circular saw blade, 165mm long
HX-GI/184a	Spare blade


Scalpel Handle

127mm long



Code	Description
HX-GI/185	No 3
HX-GI/185a	No 4
HX-GI/185b	No 4 L
HX-GI/185c	No 7

Mayo Instrument Pin



Code	Description
HX-GI/186	102mm long
HX-GI/186a	152mm long

Baron Scalpel Handle

130mm long



Code	Description
HX-GI/187	No 3
HX-GI/187a	No 4

Silicone Matting and Wire Sterilising Baskets


Silicone Matting

Medical grade for storage, sterilisation and transportation of instruments.

Cushion, protect, and organise delicate surgical instruments during storage, transportation, and sterilisation. Perforations allow for ventilation during the sterilisation process. The silicone mat can be cut to fit any size container.




SILICONE MATTING AVAILABLE FOR ALL BASKETS - PLEASE ASK FOR PRICE



502mm long x 242mm wide
7mm holes, 10mm between holes

Code


HX-GI/189



515mm long x 240mm wide
3mm holes, 6mm between holes

Code

HX-GI/189a



570mm long x 550mm wide
4mm holes, 13mm between holes

Code

HX-GI/189b

Wire Sterilising Baskets

DIN Basket

Wire mesh base with perforated sides and two carrying handles

Code	Description
HX-GI/190	230 x 130 x 60mm
HX-GI/190a	240 x 130 x 50mm
HX-GI/190b	240 x 130 x 60mm
HX-GI/190c	240 x 250 x 50mm
HX-GI/190d	240 x 250 x 60mm
HX-GI/190e	300 x 300 x 50mm
HX-GI/190f	350 x 350 x 60mm
HX-GI/190g	480 x 250 x 50mm
HX-GI/190h	480 x 250 x 60mm

Code	Description
HX-GI/190i	480 x 250 x 70mm
HX-GI/190j	480 x 400 x 60mm
HX-GI/190k	480 x 480 x 50mm
HX-GI/190l	480 x 480 x 60mm
HX-GI/190m	500 x 300 x 50mm
HX-GI/190n	550 x 430 x 50mm
HX-GI/190o	600 x 300 x 60mm



DIN Basket

Wire mesh basket with two carrying handles

Code	Description
HX-GI/191	240 x 250 x 50mm
HX-GI/191a	240 x 250 x 60mm
HX-GI/191b	240 x 250 x 70mm
HX-GI/191c	300 x 300 x 60mm
HX-GI/191d	480 x 250 x 50mm

Code	Description
HX-GI/191e	480 x 250 x 60mm
HX-GI/191f	480 x 250 x 70mm
HX-GI/191g	480 x 250 x 100mm
HX-GI/191h	480 x 480 x 50mm
HX-GI/191i	480 x 480 x 70mm

DIN Basket Lids

For perforated and wire mesh baskets

Code	Description
HX-GI/192	240 x 250mm
HX-GI/192a	300 x 300mm
HX-GI/192b	480 x 250mm
HX-GI/192c	480 x 480mm

Did you know...

Organising and fixing instruments in place can aid the cleaning process.


Help keep your instruments in good condition for longer with our baskets and trays. Available in custom sizes, to suit your requirements perfectly.

info@hypexsurgical.com

Wire Sterilising Baskets

ENT Basket

Wire mesh base with perforated sides and two carrying handles, supplied with lid



Code	Description
HX-GI/193	480 x 250 x 160mm
HX-GI/193a	Fixing clamps

Arthroscope Basket

Wire mesh base with two fixing clamps



Code	Description
HX-GI/195	290 x 80 x 52mm

Endoscope Basket

With small lens basket

Code	Description
HX-GI/196	6400 x 150 x 77mm

Fine Mesh Basket



With hinged take off / lockable lid

Code	Description
SC-GI/CO041	120 x 155 x 40mm
SC-GI/CO044	120 x 155 x 50mm
SC-GI/CO042	235 x 155 x 40mm
SC-GI/CO045	235 x 155 x 50mm
SC-GI/CO043	375 x 260 x 40mm


With separate lid

Code	Description
SC-GI/CO038	120 x 155 x 40mm
SC-GI/CO039	235 x 155 x 40mm
SC-GI/CO040	375 x 260 x 40mm

Wire mesh base with two fixing clamps

Code	Description
HX-GI/197	S - 460 x 80 x 52mm
HX-GI/197a	D - 480 x 80 x 52mm
HX-GI/197b	D - 680 x 80 x 52mm

Sieve



Code	Description
HX-GI/198	203mm x 50mm ø




WE ALSO MAKE CUSTOM BASKETS AND TRAYS,
CONTACT US TO DISCUSS YOUR REQUIREMENTS

Hollow Ware

Kidney Dish


Stainless steel



Code	Description
HX-GI/199	152mm
HX-GI/199a	203mm
HX-GI/199b	254mm
HX-GI/199c	305mm

Lotion Bowl


Without lid, stainless steel



Code	Description
HX-GI/200	102mm
HX-GI/200a	152mm
HX-GI/200b	203mm
HX-GI/200c	254mm
HX-GI/200d	305mm

Instrument Tray

Stainless steel



Code	Description
HX-GI/201	254 x 203 x 50mm
HX-GI/201a	Lid only for 254 x 203 x 50mm
HX-GI/201b	305 x 254 x 50mm
HX-GI/201c	Lid only for 305 x 254 x 50mm
HX-GI/201d	356 x 305 x 50mm
HX-GI/201e	Lid only for 356 x 305 x 50mm

Cheattle Jar



Code	Description
HX-GI/202	Stainless steel, 203mm x 50mm

Measuring Jug




Code	Description
HX-GI/203	1 litre capacity

Petri Dish



Code	Description
HX-GI/204	100mm, stainless steel, with lid

Gallipot



Code	Description
HX-GI/205	Stainless steel, 63mm

PRECISION • INNOVATION • TRUST

Crafting world-class surgical and dental instruments engineered for accuracy, durability, and performance.

About Hypex Surgical

Based in Sialkot, Pakistan – the heart of surgical manufacturing world
Hypex Surgical specializes in surgical and dental instrument. Every instrument is crafted with precision engineering, ISO certified processes, and strict quality assurance to meet the global standards of healthcare professional worldwide.

Our Commitment:

- ✓ ISO 13485 & CE Certified Manufacturing
- ✓ Sustainable & Ethical Production Practices
- ✓ 100% Traceability and Quality Guarantee

Scan to Explore
Our Product Range



E. info@hypexsurgical.com

T. +92 332 0640444

www.hypexsurgical.com

Proudly Crafted in Sialkot, Pakistan –
Serving Healthcare Professionals Worldwide

FOLLOW US

