



HYPEX SURGICAL
I N S T R U M E N T S

UROLOGY INSTRUMENTS CATALOGUE

2026 EDITION

EST. IN SIALKOT, PAKISTAN

**LEADING WORLDWIDE MANUFACTURER
OF SURGICAL & DENTAL INSTRUMENTS**



ISO Aligned



Quality Certified



Batch Traceability



Global Standards

About Us

Hypex Surgical is a Pakistan-based manufacturer of high-quality surgical and dental instruments, built on precision, craftsmanship, and trust. With years of expertise in the medical instruments industry, we take pride in supplying reliable tools that empower healthcare professionals around the world.

Our products are crafted using German-grade stainless steel and advanced production technology, ensuring accuracy, durability, and ergonomic design for every procedure. From general surgery to specialized disciplines, Hypex Surgical instruments reflect our ongoing commitment to excellence and innovation.

We proudly serve global distributors, hospitals, and medical professionals, ensuring timely delivery and unmatched quality standards — because for us, every instrument represents a promise of life-saving precision.

Our **Mission** is to design, manufacture, and deliver world-class surgical instruments that meet the evolving needs of the healthcare industry — blending innovation, craftsmanship, and integrity to support surgeons and improve patient outcomes worldwide.

Hypex Surgical **Vision** is to be recognised globally as a trusted partner in surgical excellence, setting new benchmarks in quality, precision, and reliability. Hypex Surgical aims to continuously evolve with technology, sustainability, and customer satisfaction at its core.

How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work
- 2 Email info@hypexsurgical.com or call +92 332 0640444 with the product codes and we will provide you with a quote
- 3 Upon approval of the quote we will process your order

How To Pay

We accept Visa, Mastercard, Amex cards, paypal and bank payments.



NEED A BESPOKE INSTRUMENT?

CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

Our Services

At Hypex Surgical, we take pride in being a trusted manufacturer and global supplier of high-quality surgical and dental instruments. Each instrument we produce is engineered with precision, built to perform, and designed to last — ensuring reliability in every procedure.

MANUFACTURING EXCELLENCE

We specialize in the design and production of premium-grade surgical and dental instruments using German stainless steel and cutting-edge technology. Our manufacturing process ensures precision, consistency, and durability in every product.

CUSTOM INSTRUMENT DEVELOPMENT

We collaborate closely with surgeons, hospitals, and distributors to develop customized instruments that meet unique procedural needs. From prototyping to final production, every detail is tailored to our client's specifications.

QUALITY CONTROL & ASSURANCE

Every instrument passes through a multi-stage inspection process, including dimensional verification, hardness testing, and surface finishing checks — ensuring compliance with international standards such as ISO 13485 and CE.

GLOBAL EXPORT SERVICES

We proudly serve clients in over 40 countries, providing seamless global logistics and documentation support to ensure smooth and timely deliveries.

PRIVATE LABEL & OEM SERVICES

We provide Private Label and OEM manufacturing options for distributors and brands worldwide, enabling them to market premium instruments under their own identity — supported by our unmatched craftsmanship and reliability.

POLISHING & FINISHING OPTIONS

Choose from a wide range of finishes including matte, mirror, satin, sandblast, and gold-plated, providing both functionality and aesthetic appeal to match your market's demand.

TIMELY DELIVERY & SUPPORT


With a streamlined supply chain and experienced export team, we ensure on-time delivery, transparent communication, and dedicated after-sales support to build long-term trust.

INNOVATION & R&D

Our R&D team continuously explores modern surgical technologies and ergonomic designs to enhance performance, reduce fatigue, and align with the latest trends in global healthcare.

Urology


Retractors and Forceps	04
Urethral Sounds and Bougies	05
Urethral Bougies	06
Miscellaneous	07



Retractors and Forceps


Walton Bladder Retractor

Code	Description
HX-UR/001	Parallel action, with convex lateral blades - 25mm x 70mm and curved central blade - 70mm x 215mm



Thompson Walker Bladder Retractor

Code	Description
HX-UR/002	One posterior blade 95mm x 30mm, one anterior wall blade 175mm x 25mm tapering to 20mm, two pairs of lateral swivel blades, shallow: 55mm x 22mm, deep: 75mm x 25mm




Millin Retropubic Retractor

Self retaining

Code
HX-UR/003


Millin retractor blades

Code	Description
HX-UR/004a	Short 25 x 55mm
HX-UR/004b	Long 25 x 80mm
HX-UR/004c	Centre 115 x 445mm



Millin Vulsellum Forceps

Screw joint for capsule, 225mm long



Code	Description
HX-UR/005	8-9 teeth, jaws 10mm wide

Millin Dissecting Spring Forceps

230mm long



Code	Description
HX-UR/006	3-4 teeth

Urethral Sounds and Bougies

Clutton Urethral Sounds

French gauge tapering, 265mm long



Code	Description
HX-UR/007	Set of 12 contains SC-UR/008 to SC-UR/008k

Code	Description
HX-UR/008	6 - 10
HX-UR/008a	8 - 12
HX-UR/008b	10 - 14
HX-UR/008c	12 - 16
HX-UR/008d	14 - 18
HX-UR/008e	16 - 20

Code	Description
HX-UR/008f	18 - 22
HX-UR/008g	20 - 24
HX-UR/008h	22 - 26
HX-UR/008i	24 - 28
HX-UR/008j	26 - 30
HX-UR/008k	28 - 32

Lister Urethral Sounds

English gauge with olive ends, straight, 280mm long




Code	Description
HX-UR/009	0 - 1
HX-UR/009a	0 - 2
HX-UR/009b	0.5 - 3
HX-UR/009c	1 - 4
HX-UR/009d	2 - 5
HX-UR/009e	3 - 6

Code	Description
HX-UR/009f	4 - 7
HX-UR/009g	5 - 8
HX-UR/009h	6 - 9
HX-UR/009i	7 - 10
HX-UR/009j	8 - 11
HX-UR/009k	9 - 12

Winsbury White Urethral Bougie

French gauge, child size, 175mm long



Code	Description
HX-UR/010	1-4 to 3-6
HX-UR/010a	5-8 to 7-10
HX-UR/010b	9-12 to 11-14

Canny Ryall Meatal Urethral Bougie

English gauge, double-ended, tapering
160mm long




Code	Description
HX-UR/011	1-4 and 3-6
HX-UR/011a	5-8 and 7-10
HX-UR/011b	9-12 and 11-14
HX-UR/011c	13-16 and 15-18
HX-UR/011d	17-20 and 19-22
HX-UR/011e	21-23 and 22-24
HX-UR/011f	26-28 and 29-31

Urethral Bougies

Buckston Browne Urethral Bougie

English gauge, 295mm long


Code	Description	Code	Description
HX-UR/012	4-6	HX-UR/012g	11-13
HX-UR/012a	5-7	HX-UR/012h	12-14
HX-UR/012b	6-8	HX-UR/012i	13-15
HX-UR/012c	7-9	HX-UR/012j	14-16
HX-UR/012d	8-10	HX-UR/012k	15-17
HX-UR/012e	9-11	HX-UR/012l	16-18
HX-UR/012f	10-12		



Clutton Urethral Bougie

French gauge, 260mm long

Code	Description	Code	Description
HX-UR/013	6 - 10	HX-UR/013f	18 - 22
HX-UR/013a	8 - 12	HX-UR/013g	20 - 24
HX-UR/013b	10 - 14	HX-UR/013h	22 - 26
HX-UR/013c	12 - 16	HX-UR/013i	24 - 28
HX-UR/013d	14 - 18	HX-UR/013j	26 - 30
HX-UR/013e	16 - 20	HX-UR/013k	28 - 32



Guyon Urethral Bougie

French gauge, 275mm long

Code	Description
HX-UR/014	6-24



Lister Urethral Bougie

English gauge, 270mm long

Code	Description	Code	Description
HX-UR/015	4-7	HX-UR/015f	10-13
HX-UR/015a	5-8	HX-UR/015g	11-14
HX-UR/015b	6-9	HX-UR/015h	12-15
HX-UR/015c	7-10	HX-UR/015i	13-16
HX-UR/015d	8-11	HX-UR/015j	14-17
HX-UR/015e	9-12	HX-UR/015k	15-18



Wyndham Powell Urethral Bougie


French gauge, 170mm long

Code	Description	Code	Description	Code	Description
HX-UR/016	12	HX-UR/023	19	HX-UR/030	26
HX-UR/017	13	HX-UR/024	20	HX-UR/031	27
HX-UR/018	14	HX-UR/025	21	HX-UR/032	28
HX-UR/019	15	HX-UR/026	22	HX-UR/033	29
HX-UR/020	16	HX-UR/027	23	HX-UR/034	30
HX-UR/021	17	HX-UR/028	24	HX-UR/035	31
HX-UR/022	18	HX-UR/029	25	HX-UR/036	32




Millin Bladder Neck Spreader

Self retaining, angled jaw



Code	Description
HX-UR/037	240mm long

Sieve



Code	Description
HX-UR/038	203mm x 50mm