



HYPEX SURGICAL
I N S T R U M E N T S

PLASTIC/ORAL INSTRUMENTS CATALOGUE

2026 EDITION

EST. IN SIALKOT, PAKISTAN

**LEADING WORLDWIDE MANUFACTURER
OF SURGICAL & DENTAL INSTRUMENTS**



ISO Aligned



Quality Certified



Batch Traceability



Global Standards

About Us

Hypex Surgical is a Pakistan-based manufacturer of high-quality surgical and dental instruments, built on precision, craftsmanship, and trust. With years of expertise in the medical instruments industry, we take pride in supplying reliable tools that empower healthcare professionals around the world.

Our products are crafted using German-grade stainless steel and advanced production technology, ensuring accuracy, durability, and ergonomic design for every procedure. From general surgery to specialized disciplines, Hypex Surgical instruments reflect our ongoing commitment to excellence and innovation.

We proudly serve global distributors, hospitals, and medical professionals, ensuring timely delivery and unmatched quality standards — because for us, every instrument represents a promise of life-saving precision.

Our **Mission** is to design, manufacture, and deliver world-class surgical instruments that meet the evolving needs of the healthcare industry — blending innovation, craftsmanship, and integrity to support surgeons and improve patient outcomes worldwide.

Hypex Surgical **Vision** is to be recognised globally as a trusted partner in surgical excellence, setting new benchmarks in quality, precision, and reliability. Hypex Surgical aims to continuously evolve with technology, sustainability, and customer satisfaction at its core.

How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work
- 2 Email info@hypexsurgical.com or call +92 332 0640444 with the product codes and we will provide you with a quote
- 3 Upon approval of the quote we will process your order

How To Pay

We accept Visa, Mastercard, Amex cards, paypal and bank payments.



NEED A BESPOKE INSTRUMENT?

CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

Our Services

At Hypex Surgical, we take pride in being a trusted manufacturer and global supplier of high-quality surgical and dental instruments. Each instrument we produce is engineered with precision, built to perform, and designed to last — ensuring reliability in every procedure.

MANUFACTURING EXCELLENCE

We specialize in the design and production of premium-grade surgical and dental instruments using German stainless steel and cutting-edge technology. Our manufacturing process ensures precision, consistency, and durability in every product.

CUSTOM INSTRUMENT DEVELOPMENT

We collaborate closely with surgeons, hospitals, and distributors to develop customized instruments that meet unique procedural needs. From prototyping to final production, every detail is tailored to our client's specifications.

QUALITY CONTROL & ASSURANCE

Every instrument passes through a multi-stage inspection process, including dimensional verification, hardness testing, and surface finishing checks — ensuring compliance with international standards such as ISO 13485 and CE.

GLOBAL EXPORT SERVICES

We proudly serve clients in over 40 countries, providing seamless global logistics and documentation support to ensure smooth and timely deliveries.

PRIVATE LABEL & OEM SERVICES

We provide Private Label and OEM manufacturing options for distributors and brands worldwide, enabling them to market premium instruments under their own identity — supported by our unmatched craftsmanship and reliability.

POLISHING & FINISHING OPTIONS

Choose from a wide range of finishes including matte, mirror, satin, sandblast, and gold-plated, providing both functionality and aesthetic appeal to match your market's demand.

TIMELY DELIVERY & SUPPORT

With a streamlined supply chain and experienced export team, we ensure on-time delivery, transparent communication, and dedicated after-sales support to build long-term trust.

INNOVATION & R&D

Our R&D team continuously explores modern surgical technologies and ergonomic designs to enhance performance, reduce fatigue, and align with the latest trends in global healthcare.

A close-up photograph of a man's face, specifically the right side, with black surgical markings drawn on his skin. A hand wearing a white latex glove holds a scalpel, positioned near the man's cheek. The background is a solid orange color.


Plastic/Oral

Scissors	04
Needle Holders	05
Forceps	06
Raspatories and Elevators	07
Raspatories	08
Spring Forceps, Hooks and Picks	10
Retractors	11
Osteotomes, Chisels and Gouges	13
Miscellaneous	14

Scissors

Kilner Scissors


Rounded blades, 115mm long



Code	Description
HX-PL/001a	Straight, sharp points
HX-PL/001b	Curved, blunt points

Aufricht Dissecting Scissors


Heavy pattern, triangular blades, 140mm long



Code	Description
HX-PL/002a	Straight
HX-PL/002b	Curved

Aufricht Dissecting Scissors


Light pattern, 150mm long



Code	Description
HX-PL/003a	Straight
HX-PL/003b	Curved

Scissors


Curved, sharp points



Code	Description
HX-PL/004	130mm long

Scissors


Angular, for cutting cartilage, wire cutting notch



Code	Description
HX-PL/005	127mm long

McIndoe Scissors


Curved on flat, round tapering ends



Code	Description
HX-PL/006	190mm long

Stevens Tenotomy Scissors


Supercut, 114mm long



Code	Description
HX-PL/007a	Curved/blunt
HX-PL/007b	Curved/sharp
HX-PL/007c	Straight/blunt
HX-PL/007d	Straight/sharp

Cartilage Cutting Scissors


Fine points, serrated edges



Code	Description
HX-PL/008	190mm long

Gorney Supercut Scissors


Rhytidectomy scissor, tungsten carbide, one micro serrated and one razor sharp edge



Code	Description
HX-PL/009a	125mm long
HX-PL/009b	200mm long, straight
HX-PL/009c	200mm long, curved

Fickling (Cathcart) Scissors


Angled



Code	Description
HX-PL/010	180mm long

Foster Scissor/Needle Holder


Fenestrated jaw, 125mm long



Code	Description
HX-PL/011	Left

Lawrence Needle Holder

Fine jaw, 150mm long



Code	Description
HX-PL/017a	
HX-PL/017b	Tungsten carbide, with two gold bows

Derf Needle Holder


Tungsten carbide insert, with two gold bows



Code	Description
HX-PL/012a	115mm long
HX-PL/012b	125mm long

Reverdin Suture Needle


With automatic eye, 203mm long



Code	Description
HX-PL/018a	Fully curved
HX-PL/018b	Curved
HX-PL/018c	Slightly curved
HX-PL/018d	Straight

Masing Needle Holder


Tungsten carbide insert, one angled bow



Code	Description
HX-PL/013	130mm long

Castroviejo Needle Holder

Code	Description
HX-PL/014a	140mm long, with catch
HX-PL/014b	140mm long, without catch



Diamond dust jaws


Code	Description
HX-PL/015a	180mm long, without catch
HX-PL/015b	200mm long, without catch
HX-PL/015c	215mm long, without catch



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Tungsten carbide jaws

Code	Description
HX-PL/016a	180mm long
HX-PL/016b	200mm long
HX-PL/016c	215mm long




ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Forceps

Mustarde Suture Holding Forceps


Convex and concave jaw



Code	Description
HX-PL/019	133mm long

Walsham Septum Straightening Forceps


220mm long



Code	Description
HX-PL/023a	Straight
HX-PL/023b	Right
HX-PL/023c	Left

Fickling Artery Forceps


190mm long



Code	Description
HX-PL/020a	Serrated
HX-PL/020b	1-2 teeth


Bone Cutting Forceps

Straight, compound action



Code	Description
HX-PL/024	146mm long


Asch Septum Forceps



Code	Description
HX-PL/021	235mm long

Ward Bone Cutting Forceps


Slightly angled on flat



Code	Description
HX-PL/025	178mm long

Jansen-Middleton Septum Forceps

Compound action, 178mm long



Code	Description
HX-PL/022	5mm bite

Raspatories and Elevators


Cleft Palette Raspatory



Code	Description
HX-PL/026a	Curved downwards
HX-PL/026b	Curved upwards

Cleft Palette Raspatory

Curved upwards



Code	Description
HX-PL/030a	Left
HX-PL/030b	Right

McIndoe Raspatory

For cleft palette, spoon end, 170mm long



Code	Description
HX-PL/027a	Small
HX-PL/027b	Medium
HX-PL/027c	Large

Farabeuf Elevator

Light, 160mm long



Code	Description
HX-PL/031a	12mm curved
HX-PL/031b	8mm curved

Howarth Elevator



Code	Description
HX-PL/028	220mm long

Syme Periosteal Elevator


Double ended



Code	Description
HX-PL/032	170mm long

Cryer Elevator

150mm long



Code	Description
HX-PL/029a	Right
HX-PL/029b	Left

Coupland Elevator

Curved, 150mm long




Code	Description
HX-PL/033a	3mm wide
HX-PL/033b	3.5mm wide
HX-PL/033c	4mm wide

Raspatories

Aufricht Raspatory


9mm wide, 220mm long



Code	Description
HX-PL/034a	Up stroke
HX-PL/034b	Down stroke

Durogrip Raspatory

7mm wide, 205mm long




12/14 teeth

Code	Description
HX-PL/035a	12/14 teeth
HX-PL/035b	16/20 teeth
HX-PL/035c	24/30 teeth
HX-PL/035d	36/44 teeth
HX-PL/035e	52/62 teeth

Fomon Raspatory


6.5 and 7.5mm wide, 215mm long



Code	Description
HX-PL/036a	Chequered
HX-PL/036b	Knotched

Gallagher Raspatory


190mm long



Code
HX-PL/037

Joseph Raspatory


Chequered, 165mm long



Code	Description
HX-PL/038a	Fine, 7.5mm
HX-PL/038b	Coarse, 9mm

Lewis Raspatory

8mm wide, 175mm long





Code	Description
HX-PL/039a	Fine teeth
HX-PL/039b	Coarse teeth

Maltz Raspatory

180mm long

Code	Description
HX-PL/040a	Straight, 9mm
HX-PL/040b	Straight, 5mm
HX-PL/040c	Curved, 5mm




Miller Raspatory


180mm long




Code	Description
HX-PL/041a	4.3 & 6mm
HX-PL/041b	5 & 6mm
HX-PL/041c	5 & 8mm





4.3 & 6mm




5 & 6mm




5 & 8mm




Passivation Process



Stainless steel part
free of contaminant



The chromium in the stainless
steel reacts with oxygen in
ambient air



The chromium oxide layer
is reformed and protects
against corrosion

Did you know...

Stainless steel is self passivating which provides natural corrosion resistance.

Passivation is the natural process that occurs to a surgical instrument containing chromium. The chromium in the stainless steel naturally reacts with oxygen in the air, creating the very thin passive layer which protects the surface of the steel.

Spring Forceps, Hooks and Picks

Gillies Dissecting Forceps



Code	Description
HX-PL/042	1 x 2 teeth, 152mm long

McIndoe Dissecting Forceps

152mm long



Code	Description
HX-PL/043a	1-2 teeth
HX-PL/043b	Serrated

Waugh Dissecting Forceps

Fine, serrated



Code	Description
HX-PL/044a	152mm long
HX-PL/044b	178mm long
HX-PL/044c	203mm long

College Pattern Stopping and Dressing Forceps

152mm long



Code	Description
HX-PL/045	Serrated

Mustarde Dissecting and Suture Holding Forceps

Combined double hooked jaws




Code	Description
HX-PL/046	156mm long



Childe Approximating Forceps

With magazine for Michel clips



Code	Description
HX-PL/047	180mm long

Mitchell Trimmer

165mm long



Code
HX-PL/048

Dental Pick



Code	Description
HX-PL/049	160mm long

Gillies Skin Hook

190mm long



Code	Description
HX-PL/050a	Standard
HX-PL/050b	Half Size



Joseph Double Skin Hook

165mm long



Code
HX-PL/051


Kleinert Kutz Skin Hooks



Code	Description
HX-PL/052a	16cm large
HX-PL/052b	16cm medium
HX-PL/052c	16cm small

Alms Skin Retractor


Self retaining



Code	Description
HX-PL/053a	70mm long
HX-PL/053b	100mm long

Austin Kay Retractor

85mm long



Code	Description
HX-PL/058	16mm wide

Kilner Retractor


Double ended, cat's paw



Code	Description
HX-PL/054	152mm long

Bowdler Henry Retractor

170mm long



Code	Description
HX-PL/059	12mm wide

Kilner Alae Retractor


85mm long



Code	Description
HX-PL/055a	10mm
HX-PL/055b	13mm

Tongue Retractor


190mm long



Code	Description
HX-PL/060	12mm wide

Cheek Retractor

140mm long



Code	Description
HX-PL/056	40mm wide

Ragnell Davis Retractor

140mm long




Code
HX-PL/057

Retractors

Flap Retractor


180mm long



Code	Description
HX-PL/061	32mm wide

Kilner Cheek Retractor

Double ended, 167mm long



Code	Description
HX-PL/065a	Small, 19mm and 25mm wide
HX-PL/065b	Large, 29mm and 35mm wide

Kilner Nasal Retractor


152mm long



Code	Description
HX-PL/062	7mm wide

Laster Cheek Retractor


110mm long



Code	Description
HX-PL/066	35mm wide

Hajek Retractor


Upper lip, 165mm long



Code	Description
HX-PL/063	20mm wide

Freece Oral Retractor

150mm long



Code	Description
HX-PL/064	33mm wide

Osteotomes, Chisels and Gouges

Ward (East Grinstead) Osteotome

Square handle, 178mm long



Code	Description
HX-PL/067a	3mm wide
HX-PL/067b	5mm wide

Silver Nasal Chisel

182mm long



Code	Description
HX-PL/071a	Straight
HX-PL/071b	Curved to right
HX-PL/071c	Curved to left

Gillies (French Pattern) Osteotome

140mm long



Code	Description
HX-PL/068a	3mm wide
HX-PL/068b	5mm wide
HX-PL/068c	7mm wide
HX-PL/068d	9mm wide
HX-PL/068e	11mm wide
HX-PL/068f	16mm wide

Pterigoideus (Maxilla) Chisel



Code	Description
HX-PL/072	255mm long

Jenkins Chisel

140mm long



Code	Description
HX-PL/069a	2mm wide
HX-PL/069b	4mm wide
HX-PL/069c	6mm wide
HX-PL/069d	8mm wide
HX-PL/069e	10mm wide

Tilley Septum Gouge

V shaped, straight



Code	Description
HX-PL/073	127mm long

McIndoe Nasal Chisel


152mm long



Code	Description
HX-PL/070a	5mm wide
HX-PL/070b	7mm wide
HX-PL/070c	9mm wide
HX-PL/070d	11mm wide

Fickling Dental Mallet


191mm long



Code	Description
HX-PL/074	25mm ø, tapering head

Miscellaneous

Silver Skin Graft Knife



Code	Description
HX-PL/075	190mm long

Eckhoff Mapping Pen



Code	Description
HX-PL/076a	145mm long
HX-PL/076b	Nibs packs of 10


Sommerland Surgical Marking Pen

Adjustable four-pointed tip




Code	Description
HX-PL/077	152mm long

Doyen Mouth Gag



Code	Description
HX-PL/080a	Adult
HX-PL/080b	Child


Fergusson Ackland Mouth Gag



Code	Description
HX-PL/081a	Adult, 210mm
HX-PL/081b	Child, 155mm


Gillies Suction Tube

Control valve, 180mm long




Code	Description
HX-PL/078a	2.5mm
HX-PL/078b	3mm
HX-PL/078c	3.5mm

Dental Cartridge Syringe



Code	Description
HX-PL/082	115mm long


Watson Skin Graft Knife



Code	Description
HX-PL/079a	Right-handed
HX-PL/079b	Left-handed

Freeman Areola Nipple Sizer 1.5mm

Code	Description
HX-PL/083a	25mm ø
HX-PL/083b	40mm ø
HX-PL/083c	45mm ø
HX-PL/083d	50mm ø
HX-PL/083e	55mm ø




Unna Comedo Expressor Double Ended Curette

2mm / 1mm holes



Code	Description
HX-PL/086	135mm long

Kelsey Fry Bone Awl 148mm long



Code	Description
HX-PL/084a	Straight
HX-PL/084b	Curved

Lang Steel Scoop

Double ended



Code	Description
HX-PL/087	7mm x 4mm, 140mm long


McHardy Curette

136mm long



Code	Description
HX-PL/088	1mm, 1.5mm, and 2mm ø

Joseph Nasal Saw Blade 35mm long



MacEwen Curette

Combined with seeker, 165mm long



Code	Description
HX-PL/089a	Oval
HX-PL/089b	Round

Horizontal shaft

Code	Description
HX-PL/085a	Cranked right
HX-PL/085b	Cranked left

Vertical shaft

Code	Description
HX-PL/085c	Cranked right
HX-PL/085d	Cranked left

Nero Smooth Cut Scisso

A ceramic-coated range of scissors with unbeatable durability, life cycle cost, anti-glare, cleaning efficacy, and the smoothest ever performance.



Contact us for more details on the Nero range

info@hypexsurgical.com

PRECISION • INNOVATION • TRUST

Crafting world-class surgical and dental instruments engineered for accuracy, durability, and performance.

About Hypex Surgical

Based in Sialkot, Pakistan – the heart of surgical manufacturing world
Hypex Surgical specializes in surgical and dental instrument. Every instrument is crafted with precision engineering, ISO certified processes, and strict quality assurance to meet the global standards of healthcare professional worldwide.

Our Commitment:

- ✓ ISO 13485 & CE Certified Manufacturing
- ✓ Sustainable & Ethical Production Practices
- ✓ 100% Traceability and Quality Guarantee

Scan to Explore
Our Product Range



E. info@hypexsurgical.com

T. +92 332 0640444

www.hypexsurgical.com

Proudly Crafted in Sialkot, Pakistan –
Serving Healthcare Professionals Worldwide

FOLLOW US

