



HYPEX SURGICAL
I N S T R U M E N T S

PROCTOLOGY INSTRUMENTS CATALOGUE

2026 EDITION

EST. IN SIALKOT, PAKISTAN

**LEADING WORLDWIDE MANUFACTURER
OF SURGICAL & DENTAL INSTRUMENTS**



ISO Aligned



Quality Certified



Batch Traceability



Global Standards

About Us

Hypex Surgical is a Pakistan-based manufacturer of high-quality surgical and dental instruments, built on precision, craftsmanship, and trust. With years of expertise in the medical instruments industry, we take pride in supplying reliable tools that empower healthcare professionals around the world.

Our products are crafted using German-grade stainless steel and advanced production technology, ensuring accuracy, durability, and ergonomic design for every procedure. From general surgery to specialized disciplines, Hypex Surgical instruments reflect our ongoing commitment to excellence and innovation.

We proudly serve global distributors, hospitals, and medical professionals, ensuring timely delivery and unmatched quality standards — because for us, every instrument represents a promise of life-saving precision.

Our **Mission** is to design, manufacture, and deliver world-class surgical instruments that meet the evolving needs of the healthcare industry — blending innovation, craftsmanship, and integrity to support surgeons and improve patient outcomes worldwide.

Hypex Surgical **Vision** is to be recognised globally as a trusted partner in surgical excellence, setting new benchmarks in quality, precision, and reliability. Hypex Surgical aims to continuously evolve with technology, sustainability, and customer satisfaction at its core.

How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work
- 2 Email info@hypexsurgical.com or call +92 332 0640444 with the product codes and we will provide you with a quote
- 3 Upon approval of the quote we will process your order

How To Pay

We accept Visa, Mastercard, Amex cards, paypal and bank payments.



NEED A BESPOKE INSTRUMENT?

CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

Our Services

At Hypex Surgical, we take pride in being a trusted manufacturer and global supplier of high-quality surgical and dental instruments. Each instrument we produce is engineered with precision, built to perform, and designed to last — ensuring reliability in every procedure.

MANUFACTURING EXCELLENCE

We specialize in the design and production of premium-grade surgical and dental instruments using German stainless steel and cutting-edge technology. Our manufacturing process ensures precision, consistency, and durability in every product.

CUSTOM INSTRUMENT DEVELOPMENT

We collaborate closely with surgeons, hospitals, and distributors to develop customized instruments that meet unique procedural needs. From prototyping to final production, every detail is tailored to our client's specifications.

QUALITY CONTROL & ASSURANCE

Every instrument passes through a multi-stage inspection process, including dimensional verification, hardness testing, and surface finishing checks — ensuring compliance with international standards such as ISO 13485 and CE.

GLOBAL EXPORT SERVICES

We proudly serve clients in over 40 countries, providing seamless global logistics and documentation support to ensure smooth and timely deliveries.

PRIVATE LABEL & OEM SERVICES

We provide Private Label and OEM manufacturing options for distributors and brands worldwide, enabling them to market premium instruments under their own identity — supported by our unmatched craftsmanship and reliability.

POLISHING & FINISHING OPTIONS

Choose from a wide range of finishes including matte, mirror, satin, sandblast, and gold-plated, providing both functionality and aesthetic appeal to match your market's demand.

TIMELY DELIVERY & SUPPORT


With a streamlined supply chain and experienced export team, we ensure on-time delivery, transparent communication, and dedicated after-sales support to build long-term trust.

INNOVATION & R&D

Our R&D team continuously explores modern surgical technologies and ergonomic designs to enhance performance, reduce fatigue, and align with the latest trends in global healthcare.

Proctology


Forceps	04
Forceps and Clamps	05
Scissors and Directors	06
Retractors and Proctoscopes	07
Sponge Holder	08



Forceps

Officer Biopsy Forceps


With 8mm cup jaw, to use through a sigmoidoscope



Code	Description
HX-PR/001	Shaft length 430mm

Yeoman Basket


Jawed to use through a sigmoidoscope



Code	Description
HX-PR/002a	Shaft length 250mm
HX-PR/002b	Shaft length 355mm
HX-PR/002c	Shaft length 430mm

Lloyd-Davis Sigmoidoscope


Biopsy forceps



Code	Description
HX-PR/003a	Shaft length 350mm
HX-PR/003b	Shaft length 450mm

Walton Rectal Biopsy Forceps

With laterally angled shanks



Code	Description
HX-PR/004	Shaft length 203mm

Officer Dissecting Forceps




Code	Description
HX-PR/005a	2-3 fine teeth, 150mm long
HX-PR/005b	2-3 fine teeth, 200mm long

Code	Description
HX-PR/005c	4-5 fine teeth, 200mm long
HX-PR/005d	4-5 fine teeth, 250mm long

Lloyd-Davis Rectal Occlusion Forceps


With longitudinal serrations



Code	Description
HX-PR/006	300mm long

Finch Excision Clamp


With hollowed jaws



Code	Description
HX-PR/007	75mm blade, 215mm long

Hays Colon Resection Clamp

Angled at 70°




Code	Description
HX-PR/008a	Jaw length 80mm 270mm long
HX-PR/008b	Jaw length 100mm 285mm long

Scissors and Directors

Abel Abdominal-Perineal Scissors

Straight

Code	Description
HX-PR/009	320mm long




Goligher Scissors

Straight

Code	Description
HX-PR/010a	300mm long
HX-PR/010b	330mm long

Curved


Code	Description
HX-PR/010c	300mm long
HX-PR/010d	330mm long



Lloyd-Davis Rectal Scissors


Flattened blades, blunt points
265mm long

Code	Description
HX-PR/011a	Straight
HX-PR/011b	Curved on flat



St Mark's Rectal Dilator

Code	Description
HX-PR/012a	Small, 14mm x 19mm
HX-PR/012b	Large, 16mm x 25mm



Lockhart Mummery Fistula Director

165mm long

Straight

Code	Description
HX-PR/013	With 1mm hole



90° curve

Code	Description
HX-PR/015	With 1mm hole



45° curve

Code	Description
HX-PR/014	With 1mm hole



Retrograde


Code	Description
HX-PR/016	With 1mm hole



Retractors and Proctoscopes


Lloyd-Davis Retractor

Code	Description
HX-PR/017a	Short, without tapered end, 125mm x 55mm
HX-PR/017b	Long, with tapered end, 170mm x 65mm x 50mm



Park Combined Anal Retractor


Self retaining, with 2 sets of interchangeable blades and one adjustable centre blade



Code	Description
HX-PR/018	140mm long

Park Anal Retractor


Self retaining



Code	Description
HX-PR/020	130mm, 95mm


Kelly Proctoscope

With obdurator, lumen



Code	Description
HX-PR/019a	140mm long
HX-PR/019b	200mm long
HX-PR/019c	Tapered, 22mm ø, x 100mm


Eisenhammer Speculum



Code	Description
HX-PR/021a	Blade 70mm long x 20mm wide
HX-PR/021b	Blade 95mm long x 20mm wide

Sponge Holder

Sponge Holder With spring handle



Code	Description
HX-PR/022	430mm long

PRECISION • INNOVATION • TRUST

Crafting world-class surgical and dental instruments engineered for accuracy, durability, and performance.

About Hypex Surgical

Based in Sialkot, Pakistan – the heart of surgical manufacturing world
Hypex Surgical specializes in surgical and dental instrument. Every instrument is crafted with precision engineering, ISO certified processes, and strict quality assurance to meet the global standards of healthcare professional worldwide.

Our Commitment:

- ✓ ISO 13485 & CE Certified Manufacturing
- ✓ Sustainable & Ethical Production Practices
- ✓ 100% Traceability and Quality Guarantee

Scan to Explore
Our Product Range



E. info@hypexsurgical.com

T. +92 332 0640444

www.hypexsurgical.com

Proudly Crafted in Sialkot, Pakistan –
Serving Healthcare Professionals Worldwide

FOLLOW US

