



HYPEX SURGICAL
I N S T R U M E N T S

NEURO/SPINAL INSTRUMENTS CATALOGUE

2026 EDITION

EST. IN SIALKOT, PAKISTAN

**LEADING WORLDWIDE MANUFACTURER
OF SURGICAL & DENTAL INSTRUMENTS**



ISO Aligned



Quality Certified



Batch Traceability



Global Standards

E. info@hypexsurgical.com

T. +92 332 0640444

hypexsurgical.com

About Us

Hypex Surgical is a Pakistan-based manufacturer of high-quality surgical and dental instruments, built on precision, craftsmanship, and trust. With years of expertise in the medical instruments industry, we take pride in supplying reliable tools that empower healthcare professionals around the world.

Our products are crafted using German-grade stainless steel and advanced production technology, ensuring accuracy, durability, and ergonomic design for every procedure. From general surgery to specialized disciplines, Hypex Surgical instruments reflect our ongoing commitment to excellence and innovation.

We proudly serve global distributors, hospitals, and medical professionals, ensuring timely delivery and unmatched quality standards — because for us, every instrument represents a promise of life-saving precision.

Our **Mission** is to design, manufacture, and deliver world-class surgical instruments that meet the evolving needs of the healthcare industry — blending innovation, craftsmanship, and integrity to support surgeons and improve patient outcomes worldwide.

Hypex Surgical **Vision** is to be recognised globally as a trusted partner in surgical excellence, setting new benchmarks in quality, precision, and reliability. Hypex Surgical aims to continuously evolve with technology, sustainability, and customer satisfaction at its core.

How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work
- 2 Email info@hypexsurgical.com or call +92 332 0640444 with the product codes and we will provide you with a quote
- 3 Upon approval of the quote we will process your order

How To Pay

We accept Visa, Mastercard, Amex cards, paypal and bank payments.



NEED A BESPOKE INSTRUMENT?

CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

Our Services

At Hypex Surgical, we take pride in being a trusted manufacturer and global supplier of high-quality surgical and dental instruments. Each instrument we produce is engineered with precision, built to perform, and designed to last — ensuring reliability in every procedure.

MANUFACTURING EXCELLENCE

We specialize in the design and production of premium-grade surgical and dental instruments using German stainless steel and cutting-edge technology. Our manufacturing process ensures precision, consistency, and durability in every product.

CUSTOM INSTRUMENT DEVELOPMENT

We collaborate closely with surgeons, hospitals, and distributors to develop customized instruments that meet unique procedural needs. From prototyping to final production, every detail is tailored to our client's specifications.

QUALITY CONTROL & ASSURANCE

Every instrument passes through a multi-stage inspection process, including dimensional verification, hardness testing, and surface finishing checks — ensuring compliance with international standards such as ISO 13485 and CE.

GLOBAL EXPORT SERVICES

We proudly serve clients in over 40 countries, providing seamless global logistics and documentation support to ensure smooth and timely deliveries.

PRIVATE LABEL & OEM SERVICES

We provide Private Label and OEM manufacturing options for distributors and brands worldwide, enabling them to market premium instruments under their own identity — supported by our unmatched craftsmanship and reliability.

POLISHING & FINISHING OPTIONS

Choose from a wide range of finishes including matte, mirror, satin, sandblast, and gold-plated, providing both functionality and aesthetic appeal to match your market's demand.

TIMELY DELIVERY & SUPPORT

With a streamlined supply chain and experienced export team, we ensure on-time delivery, transparent communication, and dedicated after-sales support to build long-term trust.

INNOVATION & R&D

Our R&D team continuously explores modern surgical technologies and ergonomic designs to enhance performance, reduce fatigue, and align with the latest trends in global healthcare.

Neuro/ Spinal Surgery

Forceps	04
Scissors	06
Dissectors and Probes	07
Hooks and Curettes	08
Retractors	09
Elevators	12
Suction Tubes	13
Rongeurs	14



Forceps

Adson Dissecting Forceps

Bayonet, 1-2 teeth, jaw 2.5mm at tip



Code	Description
HX-NEU/001	191mm long

Jefferson Dressing Forceps



Code	Description
HX-NEU/004	Serrated jaws, 178mm long

Adson Dissecting Forceps

Bayonet, fine, serrated jaws, jaw 2.5mm at tip



Code	Description
HX-NEU/002	178mm long

Jefferson Forceps



Code	Description
HX-NEU/005	Fine, 1-2 teeth, 178mm long

Northfield Dressing Forceps

Extra light, serrated jaws



Code	Description
HX-NEU/003	1mm wide, 203mm long

Pennybacker Dissecting Forceps



Code	Description
HX-NEU/006	3-4 teeth, 165mm long

Forceps

Silcock Iris Forceps

102mm long



Code	Description
HX-NEU/007	1-2 teeth
HX-NEU/007a	Serrated jaws

Dott Dissecting Forceps

Fine, 1-2 teeth



Code	Description
HX-NEU/010	178mm long

Cairns Dissecting Forceps


Delicate, 2-3 teeth



Code	Description
HX-NEU/008a	165mm long
HX-NEU/008b	216mm long

Dott Dressing Forceps

Angled



Code	Description
HX-NEU/011	114mm long

Cairns (Dott) Dressing Forceps


Serrated jaws



Code	Description
HX-NEU/009a	Medium, 2mm wide, 216mm long
HX-NEU/009b	Fine, 1mm wide, 178mm long

Tilley Dressing Forceps

Angled, serrated ends




Code	Description
HX-NEU/012	89mm long

Scissors

Cairns Scissors


Angled on flat



Code	Description
HX-NEU/013	178mm long

Olivecrona Scissors


Angled, with all micro serrated blade



Code	Description
HX-NEU/015	22mm blade, 150mm shaft

McKissock Scissors


Straight



Code	Description
HX-NEU/014	146mm long


Schmieden (Taylor) Scissors

Angled, probe-pointed under blade



Code	Description
HX-NEU/016	146mm long

The Original Titanium Rhoton Dissector Set




+
19 dissectors
in a purpose
designed
case

- 1 Made from lightweight titanium for balance, durability and corrosion resistance
- 2 Design features include round, deeply knurled handles which facilitate non-slip, finger-controlled micro manipulation

Order your set today
info@hypexsurgical.com

Sweet Scissors

Delicate



Code	Description
HX-NEU/017a	7mm blade, 144mm shaft
HX-NEU/017b	13mm blade, 117mm shaft

Dissectors and Probes

Cairns Guide and Director

Malleable



Code	Description
HX-NEU/018	127mm long

Dissector
200mm long



Code	Description
HX-NEU/019a	Soft
HX-NEU/019b	Soft, with hole

Pennybacker Probe Dissector

Double ended



Code	Description
HX-NEU/021	229mm long

Swedish Dissector
Double ended



Code	Description
HX-NEU/022a	Light model, 178mm long
HX-NEU/022b	Heavy model, 200mm long

Penfield Dissector



No 1

Broad curve, with 7.5mm sharp round spoon, 184mm long



No 2

7.5mm sharp round spoon with wax packer, 178mm long



No 3

Full curve and wax packer, 197mm long



No 4

Slight curve, 4mm wide, 216mm long

Code	Description
HX-NEU/020a	No 1
HX-NEU/020b	No 2
HX-NEU/020c	No 3
HX-NEU/020d	No 4

Watson Cheyne Dissector

Fine, double-ended, with probe end



Code	Description
HX-NEU/023a	127mm long
HX-NEU/023b	191mm long
HX-NEU/023c	230mm long
HX-NEU/023d	255mm long

Hooks and Curettes

Hooks

Adson Hook


178mm long



Code	Description
HX-NEU/024a	Angled
HX-NEU/024b	Straight
HX-NEU/024c	Straight, blunt

Adson Oliver Neuro Hook

191mm long



Code	Description
HX-NEU/025a	Angled downward, blunt
HX-NEU/025b	Sharp
HX-NEU/025c	Angled upward, blunt

Cairns Tissue Hook

127mm long



Code	Description
HX-NEU/026a	Half round, sharp
HX-NEU/027b	Angled, sharp

Oliver Tissue Hook

Blunt



Code	Description
HX-NEU/028	127mm long

Curettes

Falconer Scoop Curette

Double ended, angled ends



Code	Description
HX-NEU/029	191mm long

American Spinal Fusion Curettes

Blue acetal handles, 254mm long



Straight

Code	Description
HX-NEU/030a	Size 000 2mm x 4mm
HX-NEU/030b	Size 00 2.5mm x 4.5mm
HX-NEU/030c	Size 0 2.75mm x 5mm
HX-NEU/030d	Size 1 3mm x 5.5mm
HX-NEU/030e	Size 2 3.25mm x 6mm
HX-NEU/030f	Size 3 4mm x 7mm
HX-NEU/030g	Size 4 6mm x 9.5mm
HX-NEU/030h	Size 5 7mm x 14mm
HX-NEU/030i	Size 6 9mm x 17mm

Angled


Code	Description
HX-NEU/030j	Size 000 2mm x 4mm
HX-NEU/030k	Size 00 2.5mm x 4.5mm
HX-NEU/030l	Size 0 2.75mm x 5mm
HX-NEU/030m	Size 1 3mm x 5.5mm
HX-NEU/030n	Size 2 3.25mm x 6mm
HX-NEU/030o	Size 3 4mm x 7mm
HX-NEU/030p	Size 4 6mm x 9.5mm
HX-NEU/030q	Size 5 7mm x 14mm
HX-NEU/030r	Size 6 9mm x 17mm

**ALSO AVAILABLE WITH ACETAL SUREGRIP HANDLE
ADD SUFFIX BH TO CATALOGUE NUMBER**

Retractors

Mastoid Retractor


3 x 3 blunt hooks



Code	Description
HX-NEU/031	114mm long

Henderson Retractor


Self retaining, straight, 4 x 4 blunt pointed blades



Code	Description
HX-NEU/035	178mm long

Jansen (Shea) Retractor


Self retaining, blunt blades



Code	Description
HX-NEU/032	114mm long

Cairns Scalp Retractor


Curved, self retaining, swivel blades



Code	Description
HX-NEU/036	140mm long

Cone Retractor


Self retaining, hinged arms, for use in ventriculography



Code	Description
HX-NEU/033	165mm long

Schnitker Scalp Retractor


Curved, self retaining, blunt swivel blades



Code	Description
HX-NEU/037	140mm long

Jefferson Retractor


Self retaining, curved, blunt pointed blades



Code	Description
HX-NEU/034	140mm long

Cone Laminectomy Retractor

Hinged arms, with cam rack




Code	Description
HX-NEU/038	254mm long

Retractors

Cone Laminectomy Retractor

Hinged arms


Code	Description
HX-NEU/039a	Left, 165mm long
HX-NEU/039b	Right, 165mm long
HX-NEU/039c	Left, 254mm long
HX-NEU/039d	Right, 254mm long



Harvey-Jackson Laminectomy Retractor

Hinged arms, blades 38mm wide x 44mm deep


Code	Description
HX-NEU/040	305mm long



Mayo-Adson (Beckmann) Laminectomy Retractor

Hinged arms

Code	Description
HX-NEU/041	292mm long




Northfield (Cairns) Retractor

Hinged arms, one pair of blades

4 prong : 25mm wide x 16mm deep
and one pair of blades


5 prong: 32mm wide x 25mm deep

Code	Description
HX-NEU/042	292mm long



Meyerding Hemi-Laminectomy Retractor


Code	Description
HX-NEU/043a	18cm long 7" blade size 50 x 10mm
HX-NEU/043b	18cm long 7" blade size 50 x 20mm
HX-NEU/043c	18cm long 7" blade size 60 x 10mm
HX-NEU/043d	18cm long 7" blade size 60 x 20mm
HX-NEU/043e	18cm long 7" blade size 70 x 10mm
HX-NEU/043f	18cm long 7" blade size 70 x 20mm
HX-NEU/043g	18cm long blade left 90 x 10mm
HX-NEU/043h	18cm long blades left 90 x 20mm
HX-NEU/043i	18cm long blades right 90 x 10mm
HX-NEU/043j	18cm long blades right 90 x 20mm



Retractors

Cairns-Northfield Laminectomy Retractor


203mm arm length with three pairs of interchangeable swivelling blades, 76mm, 64mm, and 51mm



Code
HX-NEU/044

Hibbs Retractor


203mm long



Code	Description
HX-NEU/048a	Blades 10mm x 38mm
HX-NEU/048b	Blades 16mm x 54mm

Krayenbuhl Retractor


6mm wide



Code	Description
HX-NEU/049a	210mm long
HX-NEU/049b	240mm long

Cairns 'Rake' Scalp Retractor


Sharp, 210mm long



Code	Description
HX-NEU/045	3 prong, 16mm wide
HX-NEU/045a	4 prong, 19mm wide

Love Retractor


6mm wide



Code	Description
HX-NEU/050	210mm long

Cairns (Cushing) Scalp Retractor

Solid blade



Code	Description
HX-NEU/046	16mm wide, 210mm long

Love Retractor

203mm long




Code	Description
HX-NEU/051a	Straight, blade 5mm x 5mm
HX-NEU/051b	Shaft angled 90° at handle
HX-NEU/051c	Shaft angled 45° at handle

Cairns Strip Retractor

Malleable, 203mm long

8 sizes: 6, 8, 10, 11, 13, 16, 19, 25mm



Code	Description
HX-NEU/047a	20 SWG
HX-NEU/047b	18 SWG

O'Connell Retractor

Double ended




Code	Description
HX-NEU/052	216mm long

Elevators

Cobb Spinal Elevator

13mm wide, knurled handle



Cobb Spinal Elevator

With blue acetal handles, 279mm long


Code	Description
HX-NEU/053a	10mm wide
HX-NEU/053b	13mm wide
HX-NEU/053c	19mm wide



ALSO AVAILABLE WITH ACETAL SUREGRIP HANDLE
ADD SUFFIX BH TO CATALOGUE NUMBER

Adson Periosteal Elevator

Serrated grip, 165mm long



straight

Code	Description
HX-NEU/054	7mm wide

curved

Code	Description
HX-NEU/054a	5mm wide
HX-NEU/054b	7mm wide

McKissock Spade Periosteal Elevator

Straight



Cushing Periosteal Elevator


Slightly curved end



Code	Description
HX-NEU/055	197mm long

Pennybacker Elevator

Curved, sharp, 171mm long



Code	Description
HX-NEU/057	6mm wide

Suction Tubes

Adson Frazier Suction Tube


Angled, fully insulated, 180mm long



Code	Description
HX-NEU/058	8 F.G.

Retractor Suction Tube


With stilette, 220mm long



Code	Description
HX-NEU/059	Replacement stilette

Frazier Suction Tubes


Malleable, silver, angled, 180mm long



Code	Description
HX-NEU/060a	7 F.G.
HX-NEU/060b	8 F.G.
HX-NEU/060c	9 F.G.

Dandy Rongeur


Straight



Code	Description
HX-NEU/061	3mm bite, 178mm shaft

Falconer Rongeur

Oblique jaws




Code	Description
HX-NEU/062a	Upward cutting, 3mm bite, 110mm shaft
HX-NEU/062b	Upward cutting, 5mm bite, 110mm shaft
HX-NEU/062c	Downward cutting, 3mm bite, 85mm shaft
HX-NEU/062d	Downward cutting, 5mm bite, 110mm shaft

Rongeurs

Beck Rongeur


Straight, for intervertebral disc



Code	Description
HX-NEU/063	2mm bite, 130mm shaft

Love-Grunwald Rongeur


Cup shaped, jaws 3mm x 10mm, 178mm shaft



Code	Description
HX-NEU/066a	Straight
HX-NEU/066b	Curved upward
HX-NEU/066c	Curved downward

Cushing Rongeur


130mm shaft



Code	Description
HX-NEU/064a	2mm bite
HX-NEU/064b	3mm bite
HX-NEU/064c	4mm bite
HX-NEU/064d	5mm bite
HX-NEU/064e	6mm bite

Mackay Rongeur


140°, straight, 114mm shaft



Code	Description
HX-NEU/067	5.5mm

Cushing Rongeur

130mm shaft



150° upward cutting


Code	Description
HX-NEU/065a	2mm bite
HX-NEU/065b	3mm bite
HX-NEU/065c	4mm bite
HX-NEU/065d	5mm bite

150° downward cutting

Code	Description
HX-NEU/065e	2mm bite
HX-NEU/065f	3mm bite
HX-NEU/065g	4mm bite
HX-NEU/065h	5mm bite

Northfield Rongeur

Curved, upward cutting




Code	Description
HX-NEU/068	4mm bite, 133mm shaft

Rongeurs

Pennybacker Rongeur


120mm shaft



Code	Description
HX-NEU/069	3mm bite
HX-NEU/069a	4mm bite

Daniel Rongeur

Miniature, for neuropathology



Code	Description
HX-NEU/073	1.5mm bite, 127mm long

St Bart's Rongeur


114mm shaft



Code	Description
HX-NEU/070	3mm bite

Cairns Rongeur


Fine jaws, angled on flat, serrated jaws, curved handles



Code	Description
HX-NEU/074	4mm bite 152mm long

St Bart's Rongeur


140°, downward cutting, 114mm shaft



Code	Description
HX-NEU/071	3mm

Danish Pattern Rongeur

Angled on flat, fine jaws, curved handles, screw joint



Code	Description
HX-NEU/075	4mm bite, 203mm long

St Bart's Rongeur


140°, upward cutting, 114mm shaft



Code	Description
HX-NEU/072	3mm

Falconer Rongeur

Angled to side, light, screw joint





Code	Description
HX-NEU/076	6mm bite, 203mm long

Rongeurs

Luer-Jansen Rongeur


Light, curved on flat, compound action



Code	Description
HX-NEU/077	4mm bite, 197mm long

Olivecrona Rongeur


Angled to side, compound action, 229mm long



Code	Description
HX-NEU/078a	5mm bite
HX-NEU/078b	6mm bite

Ruskin Rongeur

Curved on flat, compound action




Code	Description
HX-NEU/079	6mm bite, 184mm long

Rongeurs

Wilms Rongeur


Angled on flat, narrow jaws



Code	Description
HX-NEU/080	4mm bite, 203mm long

Pennybacker Rongeur


Straight, fine jaws, 203mm long



Code	Description
HX-NEU/084	4mm bite

Glasgow Pattern Rongeur


Slightly curved on flat, screw joint



Code	Description
HX-NEU/081	4mm bite, 146mm long

Northfield Rongeur


Straight, 10mm bite



Code	Description
HX-NEU/085	One shallow blade, 200mm long

Heanley Rongeur


Long fine straight jaws, screw joint



Code	Description
HX-NEU/082	4mm bite, 178mm long

Stookey Rongeur


Angled to side, fine jaws



Code	Description
HX-NEU/086	3mm bite, 178mm long

Jansen Rongeur


Short curve, screw joint



Code	Description
HX-NEU/083	8mm bite, 152mm long

Trotter Rongeur

Angled to side, strong, screw joint




Code	Description
HX-NEU/087	10mm bite, 203mm long


Rongeurs

Leksell Rongeur

Curved to side, compound action




Code	Description
HX-NEU/088	Small, 6mm bite, 229mm long
HX-NEU/088a	Large, 8mm bite, 254mm long




Luer Stille Bone Rongeur

Curved, 250mm long




Code	Description
HX-NEU/089	11mm bite
HX-NEU/089a	Heavy, 229mm long



Swedish Pattern Laminectomy Shears

Double crank, compound action

Code	Description
HX-NEU/090	267mm long



Rongeurs

Kerrison Rongeurs

Cloward disc punch (not illustrated)
130° upward cutting, 190mm long




**BLACK CERAMIC COATED KERRISONS AVAILABLE
PLEASE ADD SUFFIX -BK WHEN ORDERING**

Code	Description
HX-NEU/091a	2mm wide
HX-NEU/091b	3mm wide
HX-NEU/091c	4mm wide
HX-NEU/091d	5mm wide
HX-NEU/091e	6mm wide



Kerrison punch

1mm 2mm 3mm 4mm 5mm



Non-Detachable

90° upward cutting

Code	Description
HX-NEU/092	180mm(L) x 1mm(W)
HX-NEU/093	180mm(L) x 2mm(W)
HX-NEU/094	180mm(L) x 3mm(W)
HX-NEU/095	180mm(L) x 4mm(W)
HX-NEU/096	180mm(L) x 5mm(W)
HX-NEU/097	230mm(L) x 2mm(W)
HX-NEU/098	230mm(L) x 3mm(W)
HX-NEU/099	230mm(L) x 4mm(W)
HX-NEU/100	230mm(L) x 5mm(W)

90° downward cutting

Code	Description
HX-NEU/101	180mm(L) x 1mm(W)
HX-NEU/102	180mm(L) x 2mm(W)
HX-NEU/103	180mm(L) x 3mm(W)
HX-NEU/104	180mm(L) x 4mm(W)
HX-NEU/105	180mm(L) x 5mm(W)
HX-NEU/106	230mm(L) x 2mm(W)
HX-NEU/107	230mm(L) x 3mm(W)
HX-NEU/108	230mm(L) x 4mm(W)
HX-NEU/109	230mm(L) x 5mm(W)

130° upward cutting

Code	Description
HX-NEU/128	180mm(L) x 1mm(W)
HX-NEU/129	180mm(L) x 2mm(W)
HX-NEU/130	180mm(L) x 3mm(W)
HX-NEU/131	180mm(L) x 4mm(W)
HX-NEU/132	180mm(L) x 5mm(W)
HX-NEU/133	200mm(L) x 3mm(W)
HX-NEU/134	200mm(L) x 4mm(W)
HX-NEU/135	200mm(L) x 5mm(W)
HX-NEU/136	200mm(L) x 6mm(W)
HX-NEU/137	230mm(L) x 2mm(W)
HX-NEU/138	230mm(L) x 3mm(W)
HX-NEU/139	230mm(L) x 4mm(W)
HX-NEU/140	230mm(L) x 5mm(W)

130° downward cutting

Code	Description
HX-NEU/141	180mm(L) x 1mm(W)
HX-NEU/142	180mm(L) x 2mm(W)
HX-NEU/143	180mm(L) x 3mm(W)
HX-NEU/144	180mm(L) x 4mm(W)
HX-NEU/145	180mm(L) x 5mm(W)
HX-NEU/146	230mm(L) x 2mm(W)
HX-NEU/147	230mm(L) x 3mm(W)
HX-NEU/148	230mm(L) x 4mm(W)
HX-NEU/149	230mm(L) x 5mm(W)

Detachable

90° upward cutting

Code	Description
HX-NEU/110	180(L) x 1mm(W)
HX-NEU/111	180(L) x 2mm(W)
HX-NEU/112	180(L) x 3mm(W)
HX-NEU/113	180(L) x 4mm(W)
HX-NEU/114	180(L) x 5mm(W)
HX-NEU/115	230(L) x 2mm(W)
HX-NEU/116	230(L) x 3mm(W)
HX-NEU/117	230(L) x 4mm(W)
HX-NEU/118	230(L) x 5mm(W)

90° downward cutting

Code	Description
HX-NEU/119	180(L) x 1mm(W)
HX-NEU/120	180(L) x 2mm(W)
HX-NEU/121	180(L) x 3mm(W)
HX-NEU/122	180(L) x 4mm(W)
HX-NEU/123	180(L) x 5mm(W)
HX-NEU/124	230(L) x 2mm(W)
HX-NEU/125	230(L) x 3mm(W)
HX-NEU/126	230(L) x 4mm(W)
HX-NEU/127	230(L) x 5mm(W)

130° upward cutting

Code	Description
HX-NEU/150	180(L) x 1mm(W)
HX-NEU/151	180(L) x 2mm(W)
HX-NEU/152	180(L) x 3mm(W)
HX-NEU/153	180(L) x 4mm(W)
HX-NEU/154	180(L) x 5mm(W)
HX-NEU/155	200(L) x 3mm(W)
HX-NEU/156	200(L) x 4mm(W)
HX-NEU/157	200(L) x 5mm(W)
HX-NEU/158	200(L) x 6mm(W)
HX-NEU/159	230(L) x 2mm(W)
HX-NEU/160	230(L) x 3mm(W)
HX-NEU/161	230(L) x 4mm(W)
HX-NEU/162	230(L) x 5mm(W)

130° downward cutting


Code	Description
HX-NEU/163	180(L) x 1mm(W)
HX-NEU/164	180(L) x 2mm(W)
HX-NEU/165	180(L) x 3mm(W)
HX-NEU/166	180(L) x 4mm(W)
HX-NEU/167	180(L) x 5mm(W)
HX-NEU/168	230(L) x 2mm(W)
HX-NEU/169	230(L) x 3mm(W)
HX-NEU/170	230(L) x 4mm(W)
HX-NEU/171	230(L) x 5mm(W)

Your biggest Kerrison challenges

Your biggest Kerrison challenges solved with one product:

- Productive bite every time
- Never have to use a dull Kerrison again
- Minimise the risk of dural tears
- Reduce hand and arm fatigue
- Easily opened for cleaning
- No need to open multiple trays

Contact us for more details info@hypexsurgical.com



PRECISION • INNOVATION • TRUST

Crafting world-class surgical and dental instruments engineered for accuracy, durability, and performance.

About Hypex Surgical

Based in Sialkot, Pakistan – the heart of surgical manufacturing world
Hypex Surgical specializes in surgical and dental instrument. Every instrument is crafted with precision engineering, ISO certified processes, and strict quality assurance to meet the global standards of healthcare professional worldwide.

Our Commitment:

- ✓ ISO 13485 & CE Certified Manufacturing
- ✓ Sustainable & Ethical Production Practices
- ✓ 100% Traceability and Quality Guarantee

Scan to Explore
Our Product Range



E. info@hypexsurgical.com

T. +92 332 0640444

www.hypexsurgical.com

Proudly Crafted in Sialkot, Pakistan –
Serving Healthcare Professionals Worldwide

FOLLOW US

