



HYPEX SURGICAL
INSTRUMENTS

OPHTHALMIC/ MICRO INSTRUMENTS CATALOGUE

2026 EDITION

EST. IN SIALKOT, PAKISTAN

**LEADING WORLDWIDE MANUFACTURER
OF SURGICAL & DENTAL INSTRUMENTS**



ISO Aligned



Quality Certified



Batch Traceability



Global Standards

About Us

Hypex Surgical is a Pakistan-based manufacturer of high-quality surgical and dental instruments, built on precision, craftsmanship, and trust. With years of expertise in the medical instruments industry, we take pride in supplying reliable tools that empower healthcare professionals around the world.

Our products are crafted using German-grade stainless steel and advanced production technology, ensuring accuracy, durability, and ergonomic design for every procedure. From general surgery to specialized disciplines, Hypex Surgical instruments reflect our ongoing commitment to excellence and innovation.

We proudly serve global distributors, hospitals, and medical professionals, ensuring timely delivery and unmatched quality standards — because for us, every instrument represents a promise of life-saving precision.

Our **Mission** is to design, manufacture, and deliver world-class surgical instruments that meet the evolving needs of the healthcare industry — blending innovation, craftsmanship, and integrity to support surgeons and improve patient outcomes worldwide.

Hypex Surgical **Vision** is to be recognised globally as a trusted partner in surgical excellence, setting new benchmarks in quality, precision, and reliability. Hypex Surgical aims to continuously evolve with technology, sustainability, and customer satisfaction at its core.

How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work
- 2 Email info@hypexsurgical.com or call +92 332 0640444 with the product codes and we will provide you with a quote
- 3 Upon approval of the quote we will process your order

How To Pay

We accept Visa, Mastercard, Amex cards, paypal and bank payments.



NEED A BESPOKE INSTRUMENT?

CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

Our Services

At Hypex Surgical, we take pride in being a trusted manufacturer and global supplier of high-quality surgical and dental instruments. Each instrument we produce is engineered with precision, built to perform, and designed to last — ensuring reliability in every procedure.

MANUFACTURING EXCELLENCE

We specialize in the design and production of premium-grade surgical and dental instruments using German stainless steel and cutting-edge technology. Our manufacturing process ensures precision, consistency, and durability in every product.

CUSTOM INSTRUMENT DEVELOPMENT

We collaborate closely with surgeons, hospitals, and distributors to develop customized instruments that meet unique procedural needs. From prototyping to final production, every detail is tailored to our client's specifications.

QUALITY CONTROL & ASSURANCE

Every instrument passes through a multi-stage inspection process, including dimensional verification, hardness testing, and surface finishing checks — ensuring compliance with international standards such as ISO 13485 and CE.

GLOBAL EXPORT SERVICES

We proudly serve clients in over 40 countries, providing seamless global logistics and documentation support to ensure smooth and timely deliveries.

PRIVATE LABEL & OEM SERVICES

We provide Private Label and OEM manufacturing options for distributors and brands worldwide, enabling them to market premium instruments under their own identity — supported by our unmatched craftsmanship and reliability.

POLISHING & FINISHING OPTIONS

Choose from a wide range of finishes including matte, mirror, satin, sandblast, and gold-plated, providing both functionality and aesthetic appeal to match your market's demand.

TIMELY DELIVERY & SUPPORT


With a streamlined supply chain and experienced export team, we ensure on-time delivery, transparent communication, and dedicated after-sales support to build long-term trust.

INNOVATION & R&D

Our R&D team continuously explores modern surgical technologies and ergonomic designs to enhance performance, reduce fatigue, and align with the latest trends in global healthcare.

Ophthalmic/ Micro


| | |
|------------------------------------|----|
| Forceps | 04 |
| Scissors | 06 |
| Needle Holders | 08 |
| Retractors, Speculums and Cannulas | 09 |
| Miscellaneous | 10 |



Forceps

Iris Forceps


1-2 teeth, 114mm long



| Code | Description |
|------------|-------------|
| HX-OP/001a | Straight |
| HX-OP/001b | Curved |

Iris Forceps


Serrated, 114mm long



| Code | Description |
|------------|-------------|
| HX-OP/002a | Straight |
| HX-OP/002b | Curved |

Iris Forceps


Extra fine, serrated, 94mm long



| Code | Description |
|-----------|-------------|
| HX-OP/003 | Curved |

Iris Fixation Forceps

Extra fine, 1-2 teeth, 102mm long



| Code | Description |
|-----------|-------------|
| HX-OP/004 | Curved |

Barraquer Corneal Forceps


76mm long



| Code | Description |
|-----------|-------------|
| HX-OP/005 | 1-2 teeth |

Epilation Forceps


89mm long



| Code | Description |
|-----------|-------------|
| HX-OP/006 | Oblique jaw |

Epilation Forceps


89mm long



| Code | Description |
|------------|-------------------------|
| HX-OP/007a | Straight ends, 3mm wide |
| HX-OP/007b | Round ends, 3mm wide |

Suture/Dissecting Forceps


102mm long



| Code | Description |
|-----------|-------------|
| HX-OP/008 | Extra fine |

Barraquer Fixation Forceps

102mm long




| Code | Description |
|-----------|-------------|
| HX-OP/009 | 9-10 teeth |

Forceps

Jameson Evans Cilia Forceps

120mm long




Code

HX-OP/010

Castroviejo Corneal Suture Forceps

94mm long



Code


HX-OP/014

Description

1-2 teeth

Tarsal Cyst Forceps

89mm long



Code


HX-OP/011

Description

Round ends

Arruga Intra Capsule Forceps

For removal of residual lens material during extra-capsula extraction, 120mm long



Code


HX-OP/015

Description

Curved tips

St Martin Corneal Suturing Forceps

1-2 teeth, 89mm long



Code


HX-OP/012

Description

For fine work

Desmarre Entropian Forceps

89mm long



Code

HX-OP/016

Description


Right eye

HX-OP/016a

Left eye

St Martin Micro Corneal Suturing Forceps

1-2 teeth, 89mm long



Code


HX-OP/013

Description

Fine pattern

Wilde Entropian Forceps

3 sizes, 89mm long



Code


HX-OP/017

Description

With screw action

Scissors


de Wecker Iris Scissors



| Code | Description |
|-----------|-------------|
| HX-OP/018 | 114mm long |


Castroviejo Scissors

102mm long



| Code | Description |
|-----------|-------------|
| HX-OP/019 | Curved |


Barraquer Iris Scissors



| Code | Description |
|-----------|-------------|
| HX-OP/020 | 62mm long |

Iris Scissors


102mm long



| Code | Description |
|------------|-------------|
| HX-OP/021 | Straight |
| HX-OP/021a | Curved |

Pooley Conjunctival Scissors


89mm long



| Code | Description |
|-----------|---------------------|
| HX-OP/022 | For fine dissection |

Thomas Strabismus Scissors

114mm long



| Code | Description |
|------------|-------------|
| HX-OP/023 | Straight |
| HX-OP/023a | Curved |

Strabismus Scissors


Supercut, 114mm long



| Code | Description |
|------------|-------------|
| HX-OP/024 | Straight |
| HX-OP/024a | Curved |

Iris Scissors

Supercut, 114mm long



| Code | Description |
|------------|-------------|
| HX-OP/025 | Straight |
| HX-OP/025a | Curved |

Scissors

Westcott Scissors

120mm long





| Code | Description |
|------------|-----------------|
| HX-OP/026 | Curved, blunt |
| HX-OP/026a | Curved, sharp |
| HX-OP/026b | Straight, blunt |
| HX-OP/026c | Straight, sharp |

Vannas Scissors

Scissors

| Code | Description |
|------------|--------------------------|
| HX-OP/027 | Blunt, curved 45°, 110mm |
| HX-OP/027a | Curved 85mm |
| HX-OP/027b | Straight 85mm |
| HX-OP/027c | Strongly curved 85mm |




Micro Scissors

| Code | Description |
|------------|--|
| HX-OP/028 | Delicate, straight, 85mm |
| HX-OP/028a | Delicate, curved, 85mm |
| HX-OP/028b | Delicate, strongly curved in plane, 85mm |

Needle Holders

Castroviejo Needle Holder

| Code | Description |
|------------|---------------------------|
| HX-OP/029 | 140mm long, with catch |
| HX-OP/029a | 140mm long, without catch |



Diamond dust jaws


| Code | Description |
|------------|---------------------------|
| HX-OP/030 | 180mm long, without catch |
| HX-OP/030a | 200mm long, without catch |
| HX-OP/030b | 215mm long, without catch |



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Tungsten carbide jaws

| Code | Description |
|------------|-------------|
| HX-OP/031 | 180mm long |
| HX-OP/031a | 200mm long |
| HX-OP/031b | 215mm long |




ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Silcock Needle Holder

140mm long


| Code | Description |
|------------|---------------|
| HX-OP/032 | With catch |
| HX-OP/032a | Without catch |



Micro Needle Holder

Round knurled handle, 120mm long


| Code | Description |
|------------|-------------------------|
| HX-OP/033 | Curved, without catch |
| HX-OP/033a | Curved, with catch |
| HX-OP/033b | Straight, without catch |



Micro Needle Holder

Curved, without catch


| Code | Description |
|-----------|-------------|
| HX-OP/034 | 110mm long |



Retractors, Speculums and Cannulas

Lang Eyelid Retractor

76mm long



| Code | Description |
|------------|-------------|
| HX-OP/035 | Left eye |
| HX-OP/035a | Right eye |

Desmarre Eyelid Retractor


Single, 140mm long



| Code | Description |
|------------|-------------------|
| HX-OP/036 | 12.7mm wide blade |
| HX-OP/036a | 15.9mm wide blade |
| HX-OP/036b | 19mm wide blade |

Barraquer Eye Speculum


Stainless steel wire, 44mm long



| Code | Description |
|------------|-------------|
| HX-OP/037 | Standard |
| HX-OP/037a | Small |

Speculum


Wire, 30mm wide x 45mm long



| Code | Description |
|------------|-------------|
| HX-OP/038 | Titanium |
| HX-OP/038a | Stainless |

Speculum


Wire, 25mm wide x 40mm long



| Code | Description |
|------------|-------------|
| HX-OP/039 | Titanium |
| HX-OP/039a | Stainless |

Speculum


Wire, 20mm wide x 30mm long



| code | Description |
|------------|-------------|
| HX-OP/040 | Titanium |
| HX-OP/040a | Stainless |

Clark Speculum

Wire, 35mm wide x 84mm long



| Code | Description |
|------------|----------------|
| HX-OP/041 | |
| HX-OP/041a | Infant - left |
| HX-OP/041b | Infant - right |

Lang Lachrymal Cannula

38mm long



| Code | Description |
|-----------|-----------------|
| HX-OP/042 | Cannula = 1mm ø |

Cannula


Curved, 44mm long



| Code | Description |
|-----------|-------------|
| HX-OP/043 | Flat ended |

Harrison Butler Lachrymal Cannula

38mm long




| Code | Description |
|-----------|-------------------|
| HX-OP/044 | Cannula = 0.6mm ø |

Miscellaneous

Spring Bag Clip


Serrated



| Code | Description |
|-----------|-------------|
| HX-OP/045 | 26mm wide |

Castroviejo Calliper

89mm long



| Code | Description |
|-----------|-------------|
| HX-OP/50 | Stainless |
| HX-OP/50a | Titanium |

Graefe Knife


127mm long



| Code |
|-----------|
| HX-OP/046 |

Tooke Knife


122mm long



| Code | Description |
|-----------|-------------|
| HX-OP/047 | 6 sizes |

Chalazian Sharp Spoons

122mm long



| Code | Description |
|-----------|-------------|
| HX-OP/048 | 6 sizes |

Snellen Vectis

Blunt, long loop, 140mm long



| Code | Description |
|------------|-------------|
| HX-OP/049 | Straight |
| HX-OP/049a | Angled |

PRECISION • INNOVATION • TRUST

Crafting world-class surgical and dental instruments engineered for accuracy, durability, and performance.

About Hypex Surgical

Based in Sialkot, Pakistan – the heart of surgical manufacturing world
Hypex Surgical specializes in surgical and dental instrument. Every instrument is crafted with precision engineering, ISO certified processes, and strict quality assurance to meet the global standards of healthcare professional worldwide.

Our Commitment:

- ✓ ISO 13485 & CE Certified Manufacturing
- ✓ Sustainable & Ethical Production Practices
- ✓ 100% Traceability and Quality Guarantee

Scan to Explore
Our Product Range



E. info@hypexsurgical.com

T. +92 332 0640444

www.hypexsurgical.com

Proudly Crafted in Sialkot, Pakistan –
Serving Healthcare Professionals Worldwide

FOLLOW US

